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## **Instruction for Use TPP Serological Pipette**

*Serological pipettes are design for the dispensing of aqueous liquids such as cell culture media or biological liquids.*

*The color-coding on pipette and dispenser enables easy recognition.*

*The negative graduation and the patented reservoir provide more working volume than conventional serological pipettes.*

*The specially designed dispenser with perforated peel-off tab serves as a dispenser for easy collection directly at the laboratory bench.*

### **Note:**

*Follow national regulations when handling biological material and chemical substances, use appropriate protective clothing*

*Take note of the regulations of aseptic working during the work process.*

*Do not use the mouth to pipette; always employ the pipetting-aid to fill and to dispense.*

### **Instructions:**

- *Check the packaging material and the product for intactness. Use only faultless material.*
- *Check the expiry date (EXP) on label and packaging. Only use products with valid EXP.*
- *Open the package in a sterile environment and remove a pipette for use*
- *Tear open the paper/plastic wrapping to expose the pipet's mouthpiece.*
- *Place the mouthpiece into the pipetting-aid, keeping the pipette inside its wrapping to avoid contamination. If the pipet does not fit snugly, replace the silicone pipet holder within the pipet aid nozzle.*
- *With the pipette in position, hold the aid and then remove the wrapping.*
- *Fill the pipette to the required volume.*

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#### **Advice and Recommendations:**

- *If pipetting in a sterile environment and the pipet touches any non-sterile surface, dispose of the pipet and replace it with a sterile one.*
- *Risk of contamination in case of multiple use and/or solution change of the pipette: Always replace the pipette when changing solutions or cell suspensions to avoid cross-contamination.*
- *Always use a pipetting aid (e.g.: Turbo-Fix # 94700, Peleus ball, Aspirette) and never "mouth pipette". For sterile work in cell culture, electronic pipetting aids are the means of choice. All pipetting aids have a sterile filter in the device helps to prevent contamination.*
- *Never place the pipetting aid on the workbench with a filled and/or used pipette. Liquid could flow back into the pipetting aid, contaminate and/or damaging it.*
- *Avoid the formation of aerosols during sterile work with biological agents. Let the pipette drain slowly. Never blow out the last drop out of a mechanical pipette; otherwise, the liquid film at the tip will burst. (Aerosol formation)*
- *Non-sterile solutions can be blown out to the last drop.*
- *Always hold the pipette with the tip down and avoid any contact with your hand.*
- *The serological pipette is not suitable for highly viscous liquids such as glycerol, PEG.*

### Chemical Resistance Ink

<b>+ = Good chemical resistance</b>	<b>± = Good to limited chemical resistance</b>	<b>- = Poor chemical resistance t</b>
Substance does not cause damage	Long-time and/or continuous exposure to the substance causes minor damages	Not suitable for exposure to the substance, damage may occur (deformation, discoloration, cracking, dissolution).

Solution	+ / ± / -	Solution	+ / ± / -
Acetic acid	+	Formaldehyde	+
Acetone	-	Glutaraldehyde	-
Acetonitrile	-	Glycerol	+
Ammonia solution	+	Isoamyl alcohol	+
Caustic soda	-	Isopropanol	+
Chloroform	-	Methanol	+
Dimethyl sulfoxide	+	Nitric acid	-
DMEM	+	Sulfuric acid	-
Ethanol	+	Triton X-100	+
Fetal bovine serum	+	Water	+

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**Additional:**

*Instructions for use, chemical resistance lists, and quality certificates of the respective products are available for download on the homepage [www.tpp.ch](http://www.tpp.ch).*