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Instructions for Use of TPP Vacuum Filtration «rapid»-Filtermax

The TPP filtration system «rapid»-Filtermax is suitable for sterile filtration of cell culture media, sera and aqueous solution.

The 0.22 µm polyether sulfone (PES) membrane ensures high flow rates with low protein binding and foam formation. The «rapid»-Filtermax is available as a filtration unit (filter top incl. reservoir flask) as well as a separate filter top or flask, with a DIN GL45 thread.

The filtration system is available in volume sizes 150, 250, 500 and 1000 ml. The individual system components are all identical in construction, available separately and can be combined with each other.

The unique square foot of the reservoir ensures a safe stand during sterile filtration.

The «rapid»-Filtermax is intended for single use only.

Note:

Follow national regulations when handling biological material, use appropriate protective clothing

Take note of the regulations of aseptic working during the work process.

Instructions:

- *Check the packaging material and the product for intactness. Use only faultless material.*
- *Check the expiry date (EXP) on label and packaging. Only use products with valid EXP.*
- *Open the packaging in a sterile environment and remove a complete system for use. Only open the packaging of the screw cap before use.*
- *Check the filtration system to ensure that the top filter is still firmly connected to the reservoir. Tighten it again if necessary. Take care not to overtighten the thread.*

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Advice and Recommendations:

- The graduation serves as a reference point for the filling quantity and is not an absolute value.
- Alcohol 90 % with mechanical stress (rubbing) can dissolve or blur the white colour of the lettering.

Information

- The user should provide following materials: Vacuum source; vacuum tubing; ev. pre-filters made of fibreglass for heavy-particle solutions to be pre-filtered.
- If a solution contains heavy particles to be filtered, it perform a pre-filtering using a fiberglass system.
- Solutions with a low particle density and a high filtration throughput can be attached to another reservoir for additional filtration. Please observe the rules of aseptic work.
- When filtering very dilute solutions, binding studies should be performed beforehand.
- Once the pouch is open, keep the product in a sterile area to ensure sterility.

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- *The flask is not suitable for freezing use.*
- *Do not autoclave product.*
- *Do not use the Filtermax in direct patient care applications or diagnostic procedures; it is intended for general laboratory use only.*

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To avoid possible injury due to implosions:

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- *Do not exceed 700 mm Hg differential vacuum at 25 °C*
- *Filter top: Only use sterile plastic or glass bottles that are approved for sterile filtration under vacuum. These should be free of any damage (danger of implosion). Pay attention to the max. total pressure.*
- *Wear eye-protection and perform the filtration in special safety units.*
- *Do not use hand-held vessels while filtering under vacuum.*

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Materials

Filter-membrane	Hydrophilic polyether sulfone (PES), type TPP fast flow
Flask and top	PS
Screw cap	PE
Vacuum connector	PP

Technische Daten	Technical Data	Werte / Values
Max. Gebrauchs-Temperatur [°C]	Max temp in use [°C]	max. 45
Vakuum [mm Hg] bei 25 °C	Vacuum [mm Hg] at 25 °C	700
Totvolumen [ml]	Dead volume [mL]	2.3
Vakuumanschluss [mm]	Vacuum connection [mm]	6 – 10
Wasserflussrate bei 25 °C	Water flow rate at 25 °C	22 ml/min/cm ² 0.07 mPa (0.7 bar)
Proteinbindung [µg/cm ²]	Protein binding [µg/cm ²]	< 20

		System System	Aufsatzfilter Filter Top	Flasche Flask
Produkt Nr. Abmessungen	Product No. Measurements	99150	99155	99157
Höhe [mm]	Height [mm]	103	57	63
Breite x Länge [mm]	Width x Length [mm]	93 x 93	90 x 89	93 x 93
Filtergrösse [cm ²]	Filter size [cm ²]	49	49	---
Porengrösse [µm]	Pore size [µm]	0.22	0.22	---
Volumen [ml]	Volume [mL]	150	150	150

Produkt Nr. Abmessungen	Product No. Measurements	99250	99255	99257
Höhe [mm]	Height [mm]	143	75	87
Breite x Länge [mm]	Width x Length [mm]	93 x 93	90 x 89	93 x 93
Filtergrösse [cm ²]	Filter size [cm ²]	49	49	----
Porengrösse [µm]	Pore size [µm]	0.22	0.22	---
Volumen [ml]	Volume [mL]	250	250	250

Produkt Nr. Abmessungen	Product No. Measurements	99500	99505	99507
Höhe [mm]	Height [mm]	213	111	121
Breite x Länge [mm]	Width x Length [mm]	93 x 93	90 x 89	93 x 93
Filtergrösse [cm ²]	Filter size [cm ²]	49	49	---
Porengrösse [µm]	Pore size [µm]	0.22	0.22	---
Volumen [ml]	Volume [mL]	500	500	500

Produkt Nr. Abmessungen	Product No. Measurements	99950	99955	99957
Höhe [mm]	Height [mm]	285	143	160
Breite x Länge [mm]	Width x Length [mm]	111 x 111	108 x 108	111 x 111
Filtergrösse [cm ²]	Filter size [cm ²]	69	69	---
Porengrösse [µm]	Pore size [µm]	0.22	0.22	---
Volumen [ml]	Volume [mL]	1000	1000	1000

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Additional:

Instructions for use, chemical resistance lists, and quality certificates of the respective products are available for download on the homepage www.tpp.ch