## Gebrauchsanleitung PCV-Röhrchen (packed cell volume) TPP „easy read“ Messinstrument

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## Instruction for Use TPP PCV (packed cell volume) tube TPP „easy read“ Measuring Device

*The PCV tube, together with the "easy read" measuring device, allows easy and fast measurement of biomass in cell culture, microbiology and suspension cultures of fungi and yeasts. You need a centrifuge with swing-out rotor.*

*The PCV tube is intended for single use, the "easy read" measuring device is reusable.*

### Note:

*Follow national regulations when handling biological material, use appropriate protective clothing*


*Take note of the regulations of aseptic working during the work process.*

### Instructions PCV tube

- *Check the packaging material and the product for intactness. Use only faultless material.*
- *Use only PCV tubes with a valid expiration date (EXP).*
- *Transfer a sample of 100 - 1000 µl of a mixed cell suspension into the PCV cell measuring tube.*
- *Close the tube with the cap #87008 if necessary.*
- *Centrifuge the tube for 1 minute at 2'500 x g in a swing-out rotor for 1.5 | 2 mL vessels.*
- *Read the volume of the cell pellet as described in the TPP "easy read" measuring instrument instructions.*
- *Discard the PCV tube*

### Instruction "easy read" Measuring Device

- *After centrifugation, insert the PCV tube ① into the opening provided in the "easy read" after centrifugation.*
- *Move the magnifying glass ③ until the height of the cell pellet in the capillary is equal to the cross.*
- *Read the cell pellet volume in µl with the "easy read" measuring device.*

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### Cleaning and disinfection of the "easy read"

Please note the following when choosing the cleaning and disinfecting agent:

- Choose procedures with neutral, aluminum-compatible cleaning and disinfecting agents and fully demineralized water.
- Avoid abrasive cleaning agents or scratching metal brushes, agents containing halogen or chloride: risk of corrosion.
- Thermal cleaning (autoclaving) is not recommended.
- Remove all cleaning agent residues by rinsing thoroughly.
- After cleaning/disinfection, dry the instrument completely and store it in a dry place.

### Comments and recommendations:

- The PCV method does not provide information about the cell viability of the culture.
- A new calibration line must be determined for each cell type.
- Cell pellet cannot be removed from the capillary.

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
Fig.2: A typical calibration line for a cell line (NS0) which correlates the cell count per milliliter (determined by microscopic counting) with the biomass value of the sample expressed in percent.

**Calculation PCV (%) value:**

- The measured volume is expressed in % of the volume of the transferred suspension sample. It corresponds to the total biomass of the culture.

$$PCV(\%) = \frac{\text{volume cellpellett}}{\text{volume probe}} \times 100$$

**Calculation of cell density:**

To determine the cell density from the biomass (PCV %) a calibration line must be established.

- Take a 2 ml sample from a well mixed cell suspension.
- Carry out a double determination for the PCV value (%) and for the cell density using the counting chamber (cell number/ml).
- The frequency of sampling and measurements performed depends on the growth rate (doubling time) of the cells.
- Enter the determined mean values in a table and create a diagram (see Fig. 2).
- Determine the linear equation based on the measured values (Fig. 2)

$$y = mx + b$$

m = Increase of the straight line  
 y = PCV [%]  
 x = Cell density [cell number/mL]  
 b = ordinate section

- The cell density of a sample [cell number/ml] can be calculated from the determined straight line equation:

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#### Technical Data:

#### Materials

Measuring device aluminum  
 Magnifying lens glass

Abmessungen	Measurements	87010
Volumengraduierung [µl]	Volume graduation [µl]	0.0 – 5.0
Länge x Breite x Höhe [mm]	Length x width x height [mm]	320 x 15 x 13.5
Lupe: Länge x Breite x Höhe [mm]	Magnifying lens: Length x width x height [mm]	30 x 33 x 38

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#### Additional:

Instructions for use, chemical resistance lists, and quality certificates of the respective products are available for download on the homepage [www.tpp.ch](http://www.tpp.ch)

#### Literatur:

1. Stettler, M., Jaccard, N., Hacker, D., De Jesus, M. Wurm, F. M., Jordan, M. (2006): New disposable tubes for rapid and precise biomass assessment for suspension cultures of mammalian cells. *Biotechnol. Bioeng.* 95 (6): 1228–1233.
2. Wurm, F., Tanner, R., Jordan, M. (2008) Eine neue Methode zur Bestimmung von Biomasse in der Zellkultur. *BioSpektrum* 05.08 508-509