

PR1MA™ β -Agarase

TBE/TAE tolerant agarase

Cut your agarose-based DNA purification time in half with PR1MA™ β -agarase!

Background

β -agarase has a high tolerance to inhibitors in electrophoretic buffers such as TBE and TAE, eliminating the need to perform a buffer exchange step before digestion.

PR1MA™ β -agarase has higher thermostability than other commercially available β -agarases, with a broad range of activity between 42°C to 50°C. Higher temperature digestion reduces the chances of residual agarose gelling at the end of the reaction, resulting in higher recovery of DNA or RNA.

Bulk β -Agarase is available (contact us for details at tech@midsci.com).

Application Benefits

- Complete digestion of agarose, with no agarose fragments left after the digestion
- Obtain DNA or RNA faster - no buffer exchange needed, with direct digestion in TAE or TBE
- Broad range of activity between 42°C and 50°C

Properties

- Concentration 1000 U/mL
- Storage buffer: 50 mM Tris-HCl, pH 7.5, 50 mM KCl, 1 mM DTT, 0.1 mM EDTA, 50% glycerol
- Operating Temperature: 42-50°C
- Recommended temperature: 50°C
- Storage temperature: -20°C

**These products are intended for research use only, not for diagnostic use. The safety and efficacy of these products in diagnostic or other clinical uses has not been established.*

Protocol for isolation of DNA from a low melt agarose gel

Agarose digestion

- 1) Separate DNA by electrophoresis in a low melt agarose gel prepared in 1X TAE or 0.5X TBE and stained with ethidium bromide.
- 2) Cut out the desired DNA band with a scalpel or razor blade from the gel under long wavelength UV light or a blue light transilluminator, minimizing the agarose around the band.
- 3) Determine the weight of the slice and then cut it into smaller pieces to enhance melting.
- 4) Transfer the agarose pieces to a microcentrifuge tube.
- 5) Melt the agarose at 70°C for 10 min.
- 6) Transfer the tube to 50°C and incubate for 10 min.
- 7) Add 1-2 units of β -agarase for every 200 mg of agarose.
- 8) Mix gently and incubate at 50°C for 60 min.

DNA precipitation

- 9) Adjust the salt concentration of the solution to 0.5 M NaCl, 0.3 M NaOAc or 2.5 M NH_4OAc for DNA precipitation. Mix with 2 volumes of isopropanol, mixing the solution by gently inverting the tube multiple times. Chill the tube on ice for 15 min.
- 10) Centrifuge at 15,000 X g for 15 min.
- 11) Remove the supernatant and wash the pellet with cold 80% ethanol.
- 12) Carefully remove the ethanol from the pellet.
- 13) Dry the pellet at 50°C for 5-10 min with the cap open.
- 14) Resuspend the DNA pellet in an appropriate volume of TE buffer or low TE buffer.