

PR1MA™ Reverse Transcriptase

An RNA-dependent DNA polymerase

PR1MA™ Reverse Transcriptase is an RNA-dependent DNA polymerase that can be used for complementary DNA (cDNA) synthesis from an RNA template and is ideal for use in PCR and isothermal amplification. Thunder RT is a robust enzyme that works in a broad range of temperatures (40 - 72°C) and has RNase H activity.

Properties

- Optimal temperature: 55°C
- Heat inactivation: 80°C for 10 minutes
- 10X Isothermal buffer included. Please use supplied buffer for optimal results.
- Storage temperature: -20°C
- Can be supplied in a glycerol-free/custom buffer.

Shipping & Storage

- RT is stored at -20°C in 50% glycerol, 10 mM Tris-HCl, 100 mM KCl, 1 mM DTT, 0.1 mM EDTA, pH 7.5.
- Shipped on dry or blue ice. On arrival store at -20°C for optimum stability.

We also offer a variety of encapsulated RNA controls for common targets and would love to work with you to develop a unique control to suit your needs. Contact us at tech@midsci.com

**These products are intended for research use only, not for diagnostic use. The safety and efficacy of these products in diagnostic or other clinical uses has not been established.*

Contents

PR1MA™ Reverse Transcriptase is provided at a concentration of 150 U/uL with 10X Isothermal buffer.

Background

Reverse transcriptase is an RNA-dependent DNA polymerase that can be used for complementary DNA (cDNA) synthesis from an RNA or DNA template and is ideal for use in RT-loop-mediated isothermal amplification (RT-LAMP). RT is a robust enzyme that works in a broad range of temperatures (40 - 72°C) and has RNase H activity.

Application Notes

RT is a robust enzyme used for first-strand synthesis of complementary DNA (cDNA) from RNA or single-stranded DNA templates. It is ideally suited for RT-loop-mediated isothermal amplification (LAMP) assays.

Quality Control

- RT Unit activity: A known reverse transcriptase is used to create a standard curve with a real-time qRT-PCR assay against which the activity of this enzyme is measured.
- Purity: >95% as determined by SDS-PAGE analysis
- PR1MA™ RT is free of detectable RNase and DNase (exo- and endonuclease).
- <0.05 ng contaminating host DNA per 15 U

1X Isothermal Reaction Buffer

50 mM Tris-HCl
 75 mM KCl
 3 mM MgCl₂
 10 mM DTT
 pH 8.3 at 25°C

*These products are intended for research use only, not for diagnostic use. The safety and efficacy of these products in diagnostic or other clinical uses has not been established.

General Protocol for LAMP Reaction

LAMP Reaction Mix:

Component	Stock	Volume	Final Concentration
10x Isothermal Buffer	10X	2.5 uL	1X
MgSO ₄	100 mM	1 uL	4 mM
dNTP Mix	25 mM	1.25 uL	1.25 mM
Dye (Optional)	Variable	Variable	Variable
Primer mix	20X	1.25 uL	1X
Bts Polymerase	8 U/uL	1 uL	0.32 U/uL
Rnase Inhibitor	40 U/uL	1 uL	1.6 U/uL
RT	150 U/uL	0.1 uL	0.6 U/uL
Template RNA	Variable	Variable	Variable
Nuclease-free water		to 25 uL	
Final		25 uL	

Prepare mix in a clean, nuclease-free microcentrifuge tube and incubate at 64 - 72°C for 30 min.

1. 10X Isothermal Buffer contains 20 mM MgSO₄ (2 mM per rxn); we recommend adding 4 mM MgSO₄ (on top of the 2 mM MgSO₄ contributed by the 10X Isothermal Buffer) to start and optimize assay
2. We recommend adding 4 mM MgSO₄ (on top of the 2 mM MgSO₄ contributed by the 10X Isothermal buffer) to start and optimize your assay from there
3. Intercalating dye (such as SYTO-82, SYTO-9, EvaGreen) are recommended for real time monitoring of amplification in LAMP reactions
4. A LAMP primer mix can be prepared with all 4 or 6 (with Loop) primers. A 20X primer mix should contain: 31.2 uM FIP, 31.2 uM BIP, 4 uM F3, 4 uM B3, 15.6 uM LoopF, 15.6 uM LoopB in TE or water
5. Thunder is provided at a concentration of 150 U/uL and the enzyme can be diluted into the reaction buffer to a lower concentration based on the user's needs. We recommend 3.75 U of enzyme per reaction
6. 1 ng-1 ug total RNA

Typical cDNA Synthesis Protocol

Component	Stock	Volume	Final Concentration
Template RNA	Variable	Variable	up to 1 ug
Isothermal Buffer	10X	2.5 uL	1X
50 uM Oligo(dT) ₁₂₋₁₈ or 60 uM gene-specific primer	Variable	2.5 uL	5 or 6 uM respectively
dNTP mix	10 mM	1.25 uL	0.5 mM
RNase Inhibitor	40 U/uL	0.625 uL	1 U/uL
RT	150 U/uL	0.025 uL	0.15 U/uL
Nuclease-free water		to 25 uL	
Final		25 uL	

- Combine components described in the above table
- Incubate for 5 minutes at 25°C for primer annealing
- Incubate between 40 - 72°C for 20 minutes for cDNA synthesis
- Heat inactive at 80°C for 10 minutes

PR1MA™ Reverse Transcriptase is provided at a concentration of 150 U/uL. The enzyme can be diluted into the reaction buffer for a lower concentration. Higher or lower amounts of enzyme can be used based on needs. We recommend starting at 15 U per LAMP reaction and 3.75 U per cDNA synthesis reaction.