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The TubeSpin® bioreactor with septum combines the properties of the standard TubeSpin® bioreactor #87050 with the requirements for laboratory automation. It enables the multiple access of needles and tubes in the field of laboratory automation for the exchange of sterile medium, solutions or withdrawal of cell-based products as well as high-throughput screening, cultivation and optimization processes of prokaryotic (aerobic bacteria) and eukaryotic (yeasts, algae, HEK, CHO, Sf-9) suspension cells.

The TubeSpin® bioreactor with silicone septum is available with a filter screw cap with openings (picture 1, #1) above a PTFE filter membrane, and Barcode 128.

The TubeSpin® Bioreactor with Septum is for single use.

Note:

Follow the national regulations for handling biological material and chemical substances, use the appropriate protective clothing.

Observe the regulations of aseptic working during the work process.

Instructions:

- *Check the packaging material and the product for intactness. Use only faultless material.*
- *Check the expiry date (EXP) on label and packaging. Only use products with valid EXP.*
- *The silicone septum with cross slit is accessible for needles and tips with a diameter of <2.1 mm, fill according to normal laboratory routine.*
- *Close the filter screw cap with ventilation holes. A continuous gas exchange takes place through the integrated hydrophobic membrane [0.22 µm]. If the PTFE filter membrane gets wet, the gas exchange is reduced for a short time.*

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Remarks and recommendations:

- The graduation serves as a guide for filling quantity and is no absolute value.
- Alcohol 90 % with mechanical stress (rubbing) might dissolve or blur the white color of the print.
- For freezing TPP recommends transferring samples to the TPP cryo tubes. If TubeSpin® bioreactor storage in minus temperature range is desired, TPP recommends testing the product beforehand according to your conditions.

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Measurements

Volume grad.
Length [mm]
Diameter [mm]
Max. RCF [x g]
Form
Septa
Optimal filling volume [mL]
Shaker: Recom. amplitude [mm]
Shaker: Recom. speed [rpm]

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Additional:

Instructions for use, chemical resistance lists, and quality certificates of the respective products are available for download on the homepage www.tpp.ch.