

**SPEC SHEET**    **ABT-HC-72S**

|                      |  |
|----------------------|--|
| <b>Product Name:</b> | ABT-HC-72S   |
| <b>Features:</b>     | <ul style="list-style-type: none"> <li>• Premier Solid Door Laboratory Refrigerator</li> <li>• 72 Cu. Ft. capacity</li> <li>• Three swing solid doors, self-closing, right hinged, not reversible</li> <li>• Microprocessor Temperature Controller</li> <li>• Adjustable operating temperature range: 1°-10°C</li> <li>• Standard Digital Temperature Display and Alarm Module with battery back-up to support temperature monitoring and alarms</li> <li>• Sample simulated digital temperature display; .1° Resolution</li> <li>• Adjustable temperature display in Celsius or Fahrenheit</li> <li>• Audible and visual high/low temperature alarms</li> <li>• 2 probes; 1 in air and 1 in sample bottle</li> <li>• Sensor error alarm</li> <li>• Power failure alarm</li> <li>• Resettable 1-point Min/Max History</li> <li>• Remote alarm contacts</li> <li>• Alarm validation (audible, visual, and remote alarm)</li> <li>• Alarm mute with ring back</li> <li>• Twelve adjustable shelves (adjustable in ½" increments)</li> <li>• Forced draft air circulation</li> <li>• Keyed door lock</li> <li>• Magnetic door gasket for positive seal</li> <li>• Heavy duty door handle</li> <li>• LED interior lights are safety shielded and switch controlled</li> <li>• Swivel casters (3 ½"), Probe access port (3/8")</li> </ul> |

**SPEC SHEET**    **ABT-HC-72S**

|                      |  |
|----------------------|--|
| <b>Product Name:</b> | ABT-HC-72S   |
| <b>Features:</b>     | <ul style="list-style-type: none"> <li>• White powder-coated steel exterior and interior</li> <li>• Hydrocarbon, natural refrigerant (R290)</li> <li>• Two year parts and labor warranty, plus an additional three compressor parts warranty</li> <li>• Exterior dimensions: 81" W x 34" D x 81 ¾" H with casters; 78 ¼" without casters</li> <li>• Cycle defrost</li> <li>• 115V, 60 Hz, 7.5 Amps, 1/2 HP</li> <li>• ETL listed</li> <li>• Shipping weight: 600 lbs.</li> </ul> |