

**SPEC SHEET ABT-HC-CP-49-TS**

<b>Product Name:</b>	ABT-HC-CP-49-TS
<b>Features:</b>	<ul style="list-style-type: none"> <li>• Templog Premier Glass Door Chromatography Refrigerator</li> <li>• 49 Cu. Ft. capacity</li> <li>• Two swing glass doors, self-closing, right hinged, not reversible</li> <li>• Microprocessor temperature controller</li> <li>• Adjustable operating temperature range: 1°C - 10°C</li> <li>• Touch Screen Digital Temperature Display and Alarm Module</li> <li>• Password-protected control panel</li> <li>• Sample simulated digital temperature display; .1° Resolution, adjustable temperature display in Celsius or Fahrenheit</li> <li>• Active temperature scale bar display</li> <li>• Resettable 1-point min/max history</li> <li>• Digital 24-hour temperature chart, Temperature data logging</li> <li>• USB port for data transfer and software updates</li> <li>• Real-time operational display clock</li> <li>• 4-20mA output</li> <li>• Electronic troubleshooting guide and Electronic QR Code to access Owner's Manual &amp; technical support contact information</li> <li>• Adjustable alarm volume and Audible and visual high/low temperature alarms</li> <li>• 2 probes; 1 in air and 1 in sample bottle</li> <li>• Sensor error alarm, power failure alarm, door ajar alarm, low battery alarm</li> <li>• Operation indicator</li> <li>• Alarm validation (audible, visual, and remote alarm), Alarm mute with ring back, Remote alarm contacts</li> <li>• Eight adjustable shelves (adjustable in 1/2" increments)</li> <li>• Forced draft air circulation</li> </ul>

**SPEC SHEET ABT-HC-CP-49-TS**

<b>Product Name:</b>	ABT-HC-CP-49-TS
<b>Features:</b>	<ul style="list-style-type: none"> <li>• Keyed door locks</li> <li>• One GFCI protected duplex electrical outlet</li> <li>• Two 2" access ports with spring-loaded, self-closing gasketed metal covers</li> <li>• Magnetic door gasket for positive seal and heavy duty door handles</li> <li>• LED interior lights are safety shielded and switch controlled</li> <li>• Swivel Casters (3 1/2") Probe access port (3/4")</li> <li>• Advanced hydrocarbon detector</li> <li>• White powder-coated steel exterior and interior</li> <li>• High density urethane foam cabinet insulation</li> <li>• Hydrocarbon, natural refrigerant (R290)</li> <li>• Three year parts and labor warranty, plus an additional two year compressor parts warranty</li> <li>• Exterior dimensions: 55 1/2" W* x 34 3/4" D x 81 3/4" H with casters; 78 1/4" H without casters. *Width w/access ports closed; 63 1/2" width w/access ports open</li> <li>• Cycle defrost</li> <li>• 115V, 60 Hz, 4.5 Amps, 1/4 HP</li> <li>• ETL listed</li> <li>• Shipping weight: 538 lbs.</li> <li>• Note: This unit requires two (2) outlets</li> </ul>