

SPEC SHEET ABT-HC-CP-69-TS

Product Name:	ABT-HC-CP-69-TS
Features:	<ul style="list-style-type: none"> • Templog Premier Glass Door Chromatography Refrigerator • 69 Cu. Ft. capacity • Three sliding glass doors, self-closing (energy efficient, double pane vacuum insulated) • Microprocessor temperature controller • Adjustable operating temperature range: 1°C - 10°C • Touch Screen Digital Temperature Display and Alarm Module • Password-protected control panel • Sample simulated digital temperature display; .1° Resolution, adjustable temperature display in Celsius or Fahrenheit • Active temperature scale bar display • Resettable 1-point min/max history • Digital 24-hour temperature chart, Temperature data logging • USB port for data transfer and software updates • Real-time operational display clock • 4-20mA output • Electronic troubleshooting guide and Electronic QR Code to access Owner's Manual & technical support contact info • Adjustable alarm volume • Audible and visual high/low temperature alarms • 2 probes; 1 in air and 1 in sample bottle • Sensor error alarm; power failure alarm, door ajar alarm, low battery alarm • Operation indicator • Alarm validation (audible, visual, and remote alarm) • Alarm mute with ring back • Remote alarm contacts

SPEC SHEET ABT-HC-CP-69-TS

Product Name:	ABT-HC-CP-69-TS
Features:	<ul style="list-style-type: none"> • Twelve adjustable shelves (adjustable in ½" increments) • Forced draft air circulation • Keyed door locks • One GFCI protected duplex electrical outlet • Two 2" access ports with spring-loaded, self-closing gasketed metal covers • Magnetic door gasket for positive seal and integrated door handles • LED interior lights are safety shielded and switch controlled • Swivel casters (3 ½") probe access port (¾") • Advanced hydrocarbon detector • White powder-coated steel exterior and interior • High density urethane foam cabinet insulation • Hydrocarbon, natural refrigerant (R290) • Three year parts and labor warranty, plus an additional two year compressor parts warranty • Exterior dimensions: 82 ½" W* x 33 ¼" D x 78 ¼" H without casters; 81 ¾" H with casters. *Width w/access ports closed; 90 ½" width w/access ports open • Cycle defrost • 115V, 60 Hz, 7.5 Amps, 1/2 HP • ETL listed • Shipping weight: 738 lbs. • Note: This unit requires two (2) outlets