

SPEC SHEET **ABT-HC-CS-69**

Product Name:	ABT-HC-CS-69
Features:	<ul style="list-style-type: none"> • Standard Glass Door Chromatography Refrigerator • 69 Cu. Ft. capacity • Three sliding glass doors, self-closing (energy efficient, double pane vacuum insulated) • Microprocessor temperature controller • Digital temperature display • Audible and visual high and low temperature alarms • 2 probes; 1 in air and 1 in sample bottle • Remote alarm contacts • Adjustable operating temperature range: 1°C – 10°C • Keyed door locks • One GFCI protected duplex electrical outlet • One 2" access port with spring-loaded, self-closing gasketed metal cover • Integrated door handles and magnetic door gasket for positive seal • Twelve adjustable shelves (adjustable in ½" increments) • Forced draft air circulation • LED interior lights are safety shielded and switch controlled • Swivel Casters (3 ½") • Probe access port (3/8") • Advanced hydrocarbon detector • White powder-coated steel exterior and interior • High density urethane foam cabinet insulation • Hydrocarbon, natural refrigerant (R290) • Two year parts and labor warranty, plus an additional three year compressor parts warranty

SPEC SHEET **ABT-HC-CS-69**

Product Name:	ABT-HC-CS-69
Features:	<ul style="list-style-type: none">• Exterior dimensions: 81 3/4" W x 33 1/4" D x 78 1/4" H without casters; 81 3/4" H with casters• Cycle defrost• 115V, 60 Hz, 7.5 Amps, 1/2 HP• ETL listed• Shipping weight: 738 lbs.• Note: This unit requires two (2) outlets