

SPEC SHEET ABT-HC-FFP-17

Product Name:	ABT-HC-FFP-17
Features:	<ul style="list-style-type: none"> • Standard Natural Refrigerant Flammable Storage Freezer • 17 Cu. Ft. Capacity • One solid swing door, right hinged • Microprocessor temperature controller • Digital temperature display with sample simulated probe • Audible and visual high and low temperature alarms • Remote alarm contacts • Adjustable operating temperature range: -15°C – -25°C • Six fixed shelves • Keyed door lock • Magnetic door gasket for positive seal • Leveling legs • Probe access port (3/8") • Powder coat exterior • High density urethane foam cabinet and door insulation • Natural, Hydrocarbon Refrigerant (R290) • Two year parts and labor warranty • Exterior dimensions: 30 3/4" W x 28 5/8" D x 67 3/4"H • Note: This unit must have 4" clearance on sides and back for adequate ventilation • Manual defrost • 115V, 60 Hz, 6 Amps, 1/2 HP • ETL listed • Shipping weight: 220 lbs.