

SPEC SHEET ABT-HC-LP-12

Product Name:	ABT-HC-LP-12
Features:	<ul style="list-style-type: none"> • Premier Glass Door Laboratory Refrigerator • 12 Cu. Ft. capacity • One swing glass door, self-closing, right hinged, not reversible • Microprocessor Temperature Controller • Adjustable operating temperature range: 1°-10°C • Digital Temperature Display and Alarm Module • Sample simulated digital temperature display; .1° Resolution • Adjustable temperature display in Celsius or Fahrenheit • Audible and visual high/low temperature alarms • Sensor error alarm • Power failure alarm • Resettable 1-point Min/Max History • Remote alarm contacts • 2 probes; 1 in air and 1 in sample bottle • Alarm validation (audible, visual, and remote alarm) • Alarm mute with ring back • Three adjustable shelves (adjustable in ½" increments) • Forced draft air circulation • Keyed door lock • LED interior lights are safety shielded and switch controlled • Swivel casters (3 ½"), Probe access port (3/8") • White powder-coated steel exterior and interior • High density urethane foam cabinet insulation • Hydrocarbon, natural refrigerant (R290)

SPEC SHEET **ABT-HC-LP-12**

Product Name:	ABT-HC-LP-12
Features:	<ul style="list-style-type: none">• Two year parts and labor warranty, plus an additional three compressor parts warranty• Exterior dimensions: 25" W x 29 ¾" D x 62 ¼" H without casters; 65 ¾" H with casters• Cycle defrost• 115V, 60 Hz, 3 Amps, 1/5 HP• ETL listed• Shipping weight: 258 lbs.