

SPEC SHEET ABT-HC-LP-26

Product Name:	ABT-HC-LP-26
Features:	<ul style="list-style-type: none"> • Premier Glass Door Laboratory Refrigerator • 26 Cu. Ft. capacity • One swing glass door, self-closing, right hinged, not reversible • Microprocessor Temperature Controller • Adjustable operating temperature range: 1°-10°C • Digital Temperature Display and Alarm Module • Sample simulated digital temperature display; .1° Resolution • Adjustable temperature display in Celsius or Fahrenheit • Audible and visual high/low temperature alarms • 2 probes; 1 in air and 1 in sample bottle • Sensor error alarm • Power failure alarm • Resettable 1-point Min/Max History • Remote alarm contacts • Alarm validation (audible, visual, and remote alarm) • Alarm mute with ring back • Four adjustable shelves (adjustable in 1/2" increments) • Forced draft air circulation • Keyed door lock • Magnetic door gasket for positive seal • Heavy duty door handle • LED interior lights are safety shielded and switch controlled • Swivel casters (3 1/2"), Probe access port (3/8") • White powder-coated steel exterior and interior • High density urethane foam cabinet insulation

SPEC SHEET **ABT-HC-LP-26**

Product Name:	ABT-HC-LP-26
Features:	<ul style="list-style-type: none">• Hydrocarbon, natural refrigerant (R290)• Two year parts and labor warranty, plus an additional three compressor parts warranty• Exterior dimensions: 28 3/8" W x 36 1/2" D x 78 1/4" H without casters; 81 3/4" H with casters• Cycle defrost• 115V, 60 Hz, 3 Amps, 1/5 HP• ETL listed• Shipping weight: 355 lbs.