

SPEC SHEET **ABT-HC-LS-16**

Product Name:	ABT-HC-LS-16
Features:	<ul style="list-style-type: none"> • Standard Glass Door Laboratory Refrigerator • 16 Cu. Ft. Capacity • One swing glass door, self-closing (energy efficient, double pane vacuum insulated), right hinged, not reversible • Microprocessor temperature controller and Digital temperature display • Audible and visual high and low temperature alarms • 2 probes; 1 in air and 1 in sample bottle • Remote alarm contacts • Adjustable operating temperature: 1°C to 10°C • Three adjustable shelves (adjustable in ½" increments) • Forced draft air circulation • Keyed door lock and Heavy duty door handle • Magnetic door gasket for positive seal • Swivel casters (3 ½"); Probe access port (3/8") • LED interior lights are safety shielded and switch controlled • Powder coated steel white interior and exterior • High density urethane foam cabinet insulation • Hydrocarbon, natural refrigerant (R290) • One year parts and labor warranty, plus an additional four year compressor parts warranty • Exterior dimensions: 25" W x 29 ¾" D x 79" H with casters; 75 ½" H without casters • Cycle defrost • 115V, 60 Hz, 3 Amps, 1/5 HP • ETL listed • Shipping weight: 283 lbs.