

SPEC SHEET ABT-HC-LS-26

Product Name:	ABT-HC-LS-26
Features:	<ul style="list-style-type: none"> • Standard Glass Door Laboratory Refrigerator • 26 Cu. Ft. Capacity • One swing glass door, self-closing (energy efficient, double pane vacuum insulated), right hinged • Microprocessor temperature controller and Digital temperature display • Audible and visual high and low temperature alarms • 2 probes; 1 in air and 1 in sample bottle • Remote alarm contacts • Adjustable operating temperature: 1°C to 10°C • Four adjustable shelves (adjustable in ½" increments) • Forced draft air circulation • Keyed door lock and Heavy duty door handle • Magnetic door gasket for positive seal • Swivel casters (3 ½"); Probe access port (3/8") • LED interior lights are safety shielded and switch controlled • Powder coated steel white interior and exterior • High density urethane foam cabinet insulation • Hydrocarbon, natural refrigerant (R290) • One year parts and labor warranty, plus an additional four year compressor parts warranty • Exterior dimensions: 28 3/8" W x 36 ½" D x 81 ¾" H with casters; 78 ¼" H without casters • Cycle defrost • 115V, 60 Hz, 3 Amps, 1/5 HP • ETL listed • ENERGY STAR CERTIFIED • Shipping weight: 355 lbs.