

SPEC SHEET **ABT-HC-RFC12**

Product Name:	ABT-HC-RFC12
Features:	<ul style="list-style-type: none"> • Premier Refrigerator & Freezer Combination • 12 Cu. Ft. capacity (10.5 Ref./1.7 Freezer) • Two exterior solid swing doors (solid door refrigerator/solid door freezer), right hinged, not reversible • Microprocessor temperature controller, refrigerator and freezer • Digital temperature display, refrigerator and freezer • Adjustable operating temperature range: Refrigerator: 1°C – 10°C; Freezer: -15°C – -25°C • Audible and visual high and low temperature alarms, refrigerator and freezer • 2 probes; 1 in air and 1 in sample bottle, in both Refrigerator and Freezer • Remote alarm contacts, refrigerator and freezer • Six adjustable shelves with guard rail on back to prevent product from being pushed against back wall, Refrigerator • Two adjustable shelves with guard rail on back to prevent product from being pushed against back wall, Freezer • Bottom footed shelf to allow air circulation, Refrigerator and Freezer • Forced draft circulation Refrigerator • Keyed door locks • Magnetic door gasket for positive seal • Integrated door handles • Leveling legs • Stacking kit • Two probe access ports (3/8"); one, refrigerator/one, freezer • White powder coated exterior • HFC-free refrigerant (R600a), foam insulation and packaging

SPEC SHEET **ABT-HC-RFC12**

Product Name:	ABT-HC-RFC12
Features:	<ul style="list-style-type: none">• Two year parts and labor warranty, plus an additional three year compressor parts warranty• Exterior dimensions: 23 5/8" W x 26 1/2" D x 80 1/4" H• Manual defrost freezer / cycle defrost refrigerator• 115V, 60 Hz, 2.3 Amps refrigerator, 1 Amp freezer, 1/3 HP refrigerator, 1/8 HP freezer• UL/C-UL listed• Shipping weight: 289 lbs.• *Please note this unit requires 2 outlets.