

**SPEC SHEET**    **ABT-HC-RFC7**

<b>Product Name:</b>	ABT-HC-RFC7
<b>Features:</b>	<ul style="list-style-type: none"> <li>• Premier Refrigerator &amp; Freezer Combination</li> <li>• 7 Cu. Ft. capacity (5.2 Ref./1.7 Freezer)</li> <li>• Two exterior swing doors (glass door refrigerator/solid door freezer), right hinged, not reversible</li> <li>• Microprocessor temperature controller, refrigerator and freezer</li> <li>• Digital temperature display, refrigerator and freezer</li> <li>• Audible and visual high and low temperature alarms, refrigerator and freezer</li> <li>• 2 probes; 1 in air and 1 in sample bottle, in both Refrigerator and Freezer</li> <li>• Remote alarm contacts, refrigerator and freezer</li> <li>• Adjustable operating temperature range: Refrigerator: 1°C – 10°C; Freezer: -15°C – -25°C</li> <li>• 3 adjustable shelves refrigerator/2 adjustable shelves freezer, shelves have back rails to prevent product from being pushed against back wall</li> <li>• Footed bottom shelf refrigerator, and freezer</li> <li>• Keyed door locks</li> <li>• Magnetic door gasket for positive seal and Integrated door handles</li> <li>• Leveling legs</li> <li>• Stacking kit</li> <li>• Two probe access ports (3/8"); one, refrigerator/one, freezer</li> <li>• White powder coated exterior</li> <li>• HFC-free refrigerant (R600a), foam insulation and packaging</li> <li>• Two year parts and labor warranty, plus an additional three year compressor parts warranty</li> <li>• Exterior dimensions: 23 5/8" W x 24 1/4" D x 53 1/8" H</li> <li>• Manual defrost freezer / cycle defrost refrigerator</li> </ul>

**SPEC SHEET**    **ABT-HC-RFC7**

<b>Product Name:</b>	ABT-HC-RFC7
<b>Features:</b>	<ul style="list-style-type: none"><li>• 115V, 60 Hz, 1.3 Amps refrigerator, 1 Amps freezer, 1/6 HP refrigerator, 1/8 HP freezer</li><li>• UL/C-UL listed</li><li>• Shipping weight: 207 lbs.</li><li>• *Please note this unit requires 2 outlets.</li></ul>