

SPEC SHEET ABT-HC-RFC7SA

Product Name:	ABT-HC-RFC7SA
Features:	<ul style="list-style-type: none"> • Premier Refrigerator & Auto Defrost Freezer Combination • 5.2 cu. ft. (Refrigerator) / 1.3 cu. ft. capacity (Freezer) • Two exterior swing doors (solid door refrigerator/solid door freezer), right hinged, not reversible • Microprocessor temperature controller, refrigerator and freezer • Digital temperature display, refrigerator and freezer • Audible and visual high and low temperature alarms, refrigerator and freezer • 2 probes; 1 in air and 1 in sample bottle, in both Refrigerator and Freezer • Remote alarm contacts, refrigerator and freezer • Adjustable operating temperature range: Refrigerator: 1°C – 10°C; Freezer: -15°C – -25°C • 3 adjustable shelves (refrigerator)/2 adjustable shelves (freezer), shelves have back rails to prevent product from being pushed against back wall • Footed bottom shelf refrigerator, and freezer • Keyed door locks • Magnetic door gasket for positive seal • Integrated door handles • Leveling legs • Stacking kit • Two probe access ports (3/8"); one, refrigerator/one, freezer • White powder coated exterior • HFC-free refrigerant (R600a), foam insulation and packaging • Two year parts and labor warranty, plus an additional three year compressor parts warranty

SPEC SHEET **ABT-HC-RFC7SA**

Product Name:	ABT-HC-RFC7SA
Features:	<ul style="list-style-type: none"> • Exterior dimensions: 23 3/4" W x 24" D x 53 1/8" H • Auto defrost freezer / cycle defrost refrigerator • 115V, 60 Hz, 1.3 Amps refrigerator, 1 Amps freezer, 1/6 HP refrigerator, 1/8 HP freezer • UL/C-UL listed • Shipping weight: 215 lbs. • *Please note this unit requires 2 outlets.