

SPEC SHEET ABT-HC-RFC9G-LH

Product Name:	ABT-HC-RFC9G-LH
Features:	<ul style="list-style-type: none"> • Premier Refrigerator & Freezer Combination - LEFT HINGED • 9 Cu. Ft. Capacity (4.6 Ref./4.2 Freezer) • Two exterior swing doors (glass door refrigerator/solid door freezer), left hinged • Microprocessor temperature controller, refrigerator and freezer • Digital temperature display, refrigerator and freezer • Audible and visual high and low temperature alarms, refrigerator and freezer • 2 probes; 1 in air and 1 in sample bottle, in both Refrigerator and Freezer • Remote alarm contacts, refrigerator and freezer • Adjustable operating temperature range: Refrigerator: 1°C – 10°C; Freezer: -15°C - -25°C • 2 adjustable/1 fixed shelves refrigerator with guard rail on back to prevent product from being pushed against back wall/2 fixed shelves freezer • Keyed door locks • Magnetic door gasket for positive seal • Recessed door handle, freezer; heavy duty door handle, refrigerator • Leveling legs • Stacking kit • LED interior light, Refrigerator • Two probe access ports (3/8"); one, refrigerator/one, freezer • White powder coated exterior • High density urethane foam cabinet and door insulation • HFC-free refrigerant (R600a), foam insulation and packaging

SPEC SHEET **ABT-HC-RFC9G-LH**

Product Name:	ABT-HC-RFC9G-LH
Features:	<ul style="list-style-type: none"> • Exterior dimensions: 23 ¾" W x 25 ½" D x 66 ¾" H • Two year parts and labor warranty, plus an additional three year compressor parts warranty • Manual defrost freezer / cycle defrost refrigerator • 115V, 60 Hz, 1.74 Amps Refrigerator, 1.5 Amps Freezer, 1/8 HP Refrigerator and Freezer • UL/C-UL listed • Shipping weight: 282 lbs. • *Please note this unit requires 2 outlets.