

SPEC SHEET ABT-HC-UCBI-0404

| | |
|----------------------|---|
| Product Name: | ABT-HC-UCBI-0404 |
| Features: | <ul style="list-style-type: none"> • Premier Undercounter Refrigerator Built-In • 4.6 Cu. Ft. Capacity • One solid swing door, right hinged* • Microprocessor temperature controller • Digital temperature display • Audible and visual high and low temperature alarms • Remote alarm contacts • 2 probes; 1 in air and 1 in glycol bottle • Adjustable operating temperature range: 1°C – 10°C • Three shelves (two adjustable/one fixed) with guard rail on back to prevent product from being pushed against the back wall • Back wall evaporator cover • Forced draft circulation • Recessed door handle • Keyed door lock • Magnetic door gasket for positive seal • Leveling legs • Probe access port (3/8") • Powder coat exterior • High density urethane foam cabinet and door insulation • HFC-free refrigerant (R600a), foam insulation and packaging • Two year parts and labor warranty, plus an additional three year compressor parts only warranty • Pyxis®, Omnicell® and AcuDose RX® compatible • Exterior dimensions: 23 ¾" W x 24 ½" D x 33 3/8" H • Cycle defrost |

SPEC SHEET **ABT-HC-UCBI-0404**

| | |
|----------------------|---|
| Product Name: | ABT-HC-UCBI-0404 |
| Features: | <ul style="list-style-type: none">• 115V, 60 Hz, 1/8 HP, 1.74 Amps• Energy Consumption 1.15 kWh/day• UL/C-UL listed• ENERGY STAR CERTIFIED• Shipping weight: 137 lbs.• *This product is available left hinged, product number: ABT-HC-UCBI-0404-LH |