

**SPEC SHEET**    **ABT-HC-UCBI-0404SS-LH**

<b>Product Name:</b>	ABT-HC-UCBI-0404SS-LH
<b>Features:</b>	<ul style="list-style-type: none"> <li>• Premier Undercounter Refrigerator Built-In-LEFT HINGED</li> <li>• 4.5 Cu. Ft. Capacity</li> <li>• One solid stainless steel swing door, left hinged</li> <li>• Microprocessor temperature controller</li> <li>• Digital temperature display</li> <li>• Audible and visual high and low temperature alarms</li> <li>• 2 probes; 1 in air and 1 in sample bottle</li> <li>• Remote alarm contacts</li> <li>• Adjustable operating temperature range: 1°C – 10°C</li> <li>• Three shelves (two adjustable/one fixed) with guard rail on back to prevent product from pushing against back wall</li> <li>• Back wall evaporator cover</li> <li>• Forced draft circulation</li> <li>• Recessed door handle</li> <li>• Keyed door lock</li> <li>• Magnetic door gasket for positive seal</li> <li>• Leveling legs</li> <li>• Probe access port (3/8")</li> <li>• Stainless steel exterior</li> <li>• High density urethane foam cabinet and door insulation</li> <li>• HFC-free refrigerant (R600a), foam insulation and packaging</li> <li>• Two year parts and labor warranty, plus an additional three year compressor parts only warranty</li> <li>• Pyxis®, Omnicell® and AcuDose RX® compatible</li> <li>• Exterior dimensions: 23 ¾" W x 24 ½" D x 33 ⅜" H</li> <li>• Cycle defrost</li> </ul>

**SPEC SHEET**    **ABT-HC-UCBI-0404SS-LH**

<b>Product Name:</b>	ABT-HC-UCBI-0404SS-LH
<b>Features:</b>	<ul style="list-style-type: none"><li>• 115V, 60 Hz, 1.74 Amps, 1/8 HP</li><li>• Energy Consumption 1.15 kWh/day</li><li>• UL/C-UL listed</li><li>• ENERGY STAR CERTIFIED</li><li>• Shipping weight: 135 lbs.</li></ul>