

SPEC SHEET ABT-HC-UCBI-0420-ADA

Product Name:	ABT-HC-UCBI-0420-ADA
Features:	<ul style="list-style-type: none"> • Premier Undercounter Freezer Built-In • 4.2 Cu. Ft. Capacity • One solid swing door, right hinged* • Microprocessor temperature controller • Digital temperature display • Audible and visual high and low temperature alarms • 2 probes; 1 in air and 1 in sample bottle • Remote alarm contacts • Adjustable operating temperature range: -15°C – -25°C • Two fixed shelves • Keyed door lock • Magnetic door gasket for positive seal • Recessed door handle • Leveling legs • Probe access port (3/8") • Powder coat exterior • High density urethane foam cabinet and door insulation • HFC-free refrigerant (R600a), foam insulation and packaging • Two year parts and labor warranty, plus an additional three year compressor parts only warranty • Pyxis®, Omnicell® and AcuDose RX® compatible • Exterior dimensions: 23 ¾" W x 24 ½" D x 31 15/16" H • Manual defrost • 115V, 60 Hz, 1.5 Amps, 1/8 HP • Energy Consumption 0.97 kWh/day • ADA compliant

SPEC SHEET **ABT-HC-UCBI-0420-ADA**

Product Name:	ABT-HC-UCBI-0420-ADA
Features:	<ul style="list-style-type: none">• UL/C-UL listed• ENERGY STAR CERTIFIED• Shipping weight: 133 lbs.• *This product is available left hinged, product number: ABT-HC-UCBI-0420-ADA-LH