

**SPEC SHEET ABT-HC-UCBI-0420SS**

<b>Product Name:</b>	ABT-HC-UCBI-0420SS
<b>Features:</b>	<ul style="list-style-type: none"> <li>• Premier Undercounter Freezer Built-In</li> <li>• 4.2 Cu. Ft. Capacity</li> <li>• One stainless steel swing door, right hinged*</li> <li>• Microprocessor temperature controller</li> <li>• Digital temperature display</li> <li>• Audible and visual high and low temperature alarms and Remote alarm contacts</li> <li>• 2 probes; 1 in air and 1 in sample bottle</li> <li>• Adjustable operating temperature range: -15°C – -25°C</li> <li>• Two fixed shelves</li> <li>• Keyed door lock and Recessed door handle</li> <li>• Magnetic door gasket for positive seal</li> <li>• Leveling legs</li> <li>• Probe access port (3/8")</li> <li>• Stainless steel exterior</li> <li>• High density urethane foam cabinet and door insulation</li> <li>• HFC-free refrigerant (R600a), foam insulation and packaging</li> <li>• Two year parts and labor warranty, plus an additional three year compressor parts only warranty</li> <li>• Pyxis®, Omnicell® and AcuDose RX® compatible</li> <li>• Exterior dimensions: 23 ¾" W x 24 ½" D x 33 3/8" H</li> <li>• Manual defrost</li> <li>• 115V, 60 Hz, 1.5 Amps, 1/8 HP</li> <li>• Energy Consumption 0.97 kWh/day UL/C-UL listed</li> <li>• ENERGY STAR CERTIFIED</li> <li>• Shipping weight: 131 lbs.</li> <li>• *This product is available left hinged, product number: ABT-HC-UCBI-0420SS-LH</li> </ul>