

SPEC SHEET **ABT-HC-UCBI-0420SS-LH**

Product Name:	ABT-HC-UCBI-0420SS-LH
Features:	<ul style="list-style-type: none"> • Premier Undercounter Freezer Built-In-LEFT HINGED • 4.2 Cu. Ft. Capacity • One stainless steel swing door, left hinged • Microprocessor temperature controller and Digital temperature display • Audible and visual high and low temperature alarms • 2 probes; 1 in air and 1 in sample bottle • Remote alarm contacts • Adjustable operating temperature range: -15°C – -25°C • Two fixed shelves • Keyed door lock • Recessed door handle and Magnetic door gasket for positive seal • Leveling legs • Probe access port (3/8") • Stainless steel exterior • High density urethane foam cabinet and door insulation • HFC-free refrigerant (R600a), foam insulation and packaging • Two year parts and labor warranty, plus an additional three year compressor parts only warranty • Pyxis®, Omnicell® and AcuDose RX® compatible • Exterior dimensions: 23 ¾" W x 24 ½" D x 33 ⅜" H • Manual defrost • 115V, 60 Hz, 1.5 Amps, 1/8 HP • Energy Consumption 0.97 kWh/day UL/C-UL listed • ENERGY STAR CERTIFIED • Shipping weight: 131 lbs.