

SPEC SHEET CRT-ABT-HC-S49S

Product Name:	CRT-ABT-HC-S49S
Features:	<ul style="list-style-type: none"> • Upright Controlled Room Temperature Cabinet, Solid Door • 49 Cu. Ft. capacity • Two swing solid doors, self-closing, right hinged, not reversible • Microprocessor temperature controller • Digital temperature display • Audible and visual high and low temperature alarms • Remote alarm contacts • Adjustable operating temperature range: 68°F–77°F (20°C – 25°C) • Temperature monitor device (TMD) complies with the current CDC guidelines, with 3 years certification of calibration, “buffered” probe in the product simulated solution, min/max memory. F/C switchable, field installable, and visual & audible temp alarm • Keyed door locks • Magnetic door gasket for positive seal • 2 probes; 1 in air and 1 in sample bottle • Eight adjustable shelves • Vaccine storage power cord warning label • Forced draft air circulation • Heavy duty door handles • LED interior lights are safety shielded and switch controlled • Swivel Casters (3”) • Probe access port (3/8”) • White powder-coated steel exterior and interior • Hydrocarbon, natural refrigerant (R290)

SPEC SHEET CRT-ABT-HC-S49S

Product Name:	CRT-ABT-HC-S49S
Features:	<ul style="list-style-type: none"> • One year parts and labor warranty, plus an additional four year compressor parts warranty • Exterior dimensions: 54" W x 34 ¾" D x 81 ¾" H with casters; 78 ¼" H without casters • Cycle defrost • 115V, 60 Hz, 4.5 Amps, 1/4 HP • UL/C-UL listed • Shipping weight: 446 lbs.