

SPEC SHEET PH-ABT-HC-12G

Product Name:	PH-ABT-HC-12G
Features:	<ul style="list-style-type: none"> • Premier Pharmacy/Vaccine Glass Door Refrigerator • 12 Cu. Ft. capacity • One swing glass door, self-closing (energy efficient, double pane vacuum insulated), right hinged, not reversible • Microprocessor Temperature Controller • Adjustable operating temperature range: 36°F – 46°F (2°C – 8°C) • Temperature monitor device (TMD) complies with the current CDC guidelines, with 3 years certification of calibration, “buffered” probe in the product simulated solution, min/max memory. F/C switchable, field installable, and visual & audible temp alarm • Standard Digital Temperature Display and Alarm Module with battery back-up to support temperature monitoring and alarms • Sample simulated digital temperature display; .1° Resolution • Adjustable temperature display in Celsius or Fahrenheit • Audible and visual high/low temperature alarms • 2 probes; 1 in air and 1 in sample bottle • Sensor error alarm • Power failure alarm • Resettable 1-point Min/Max History • Remote alarm contacts • Alarm validation (audible, visual, and remote alarm) • Alarm mute with ring back • 5 total shelves, 3 adjustable plus 2 additional • Vaccine storage power cord warning label • Forced draft air circulation

SPEC SHEET PH-ABT-HC-12G

Product Name:	PH-ABT-HC-12G
Features:	<ul style="list-style-type: none"> • Keyed door lock • Magnetic door gasket for positive seal • Heavy duty door handle • Swivel Casters (3 ½") • LED interior lights are safety shielded and switch controlled • Probe access port (¾") • Powder coated steel white interior and exterior • Hydrocarbon, natural refrigerant (R290) • Pharmacy refrigerator/freezer toolkit and temperature logs • Two year parts and labor warranty, plus an additional three year compressor parts warranty • Exterior dimensions: 25" W x 29 ¾" D x 65 ¾" H with casters; 62 ¼" H without casters • Cycle defrost • 115V, 60 Hz, 3 Amps, 1/5 HP • ETL listed • Shipping weight: 258 lbs