

**SPEC SHEET PH-ABT-HC-S69G**

<b>Product Name:</b>	PH-ABT-HC-S69G
<b>Features:</b>	<ul style="list-style-type: none"> <li>• Standard Pharmacy/Vaccine Glass Door Refrigerator</li> <li>• 69 Cu. Ft. Capacity</li> <li>• Three sliding glass doors, self-closing (energy efficient, double pane vacuum insulated)</li> <li>• Microprocessor temperature controller</li> <li>• Digital temperature display</li> <li>• Audible and visual high and low temperature alarms</li> <li>• Remote alarm contacts</li> <li>• 2 probes; 1 in air and 1 in sample bottle</li> <li>• Adjustable operating temperature range: 36°F – 46°F (2°C – 8°C)</li> <li>• Refrigerator/Freezer freestanding thermometer with 3 year certificate of calibration, with min/max memory, °F/°C switchable</li> <li>• Twelve adjustable shelves (adjustable in ½" increments)</li> <li>• Vaccine storage power cord warning label</li> <li>• Forced draft air circulation</li> <li>• Keyed door lock</li> <li>• Magnetic door gasket for positive seal</li> <li>• Integrated door handles</li> <li>• Swivel Casters (3 ½")</li> <li>• LED interior lights are safety shielded and switch controlled</li> <li>• Probe access port (¾")</li> <li>• Powder coated steel white interior and exterior</li> <li>• High density urethane foam cabinet insulation</li> <li>• Hydrocarbon, natural refrigerant (R290)</li> <li>• Pharmacy refrigerator/freezer toolkit and temperature logs</li> </ul>

**SPEC SHEET PH-ABT-HC-S69G**

<p><b>Product Name:</b></p>	<p>PH-ABT-HC-S69G</p>
<p><b>Features:</b></p>	<ul style="list-style-type: none"> <li>• One year parts and labor warranty, plus an additional four year compressor parts warranty</li> <li>• Exterior dimensions: 81" W x 33 ¼" D x 81 ¾" H with casters; 78 ¼" H without casters</li> <li>• Cycle defrost</li> <li>• 115V, 60 Hz, 7.5 Amps, 1/2 HP</li> <li>• ETL listed</li> <li>• Shipping weight: 738 lbs.</li> </ul>