

**SPEC SHEET PH-ABT-HC-SSP-49**

<b>Product Name:</b>	PH-ABT-HC-SSP-49
<b>Features:</b>	<ul style="list-style-type: none"> <li>• Premier Pharmacy/Vaccine Stainless Steel Refrigerator</li> <li>• 49 Cu. Ft. Capacity</li> <li>• Two stainless steel swing doors</li> <li>• Microprocessor temperature controller</li> <li>• Digital temperature display</li> <li>• Audible and visual high and low temperature alarms</li> <li>• Remote alarm contacts</li> <li>• 2 probes; 1 in air and 1 in sample bottle</li> <li>• Adjustable operating temperature range: 36°F – 46°F (2°C- 8°C)</li> <li>• Refrigerator/Freezer freestanding thermometer with 3 year certificate of calibration, with min/max memory, °F/°C switchable</li> <li>• Six adjustable shelves</li> <li>• Vaccine storage power cord warning label</li> <li>• Forced draft air circulation</li> <li>• Keyed door locks</li> <li>• Magnetic door gasket for positive seal</li> <li>• Integrated door handles</li> <li>• Swivel casters (adds 6 1/2" to height)</li> <li>• Interior light</li> <li>• Probe access port (3/8")</li> <li>• Stainless steel interior floor</li> <li>• Stainless steel exterior</li> <li>• High density urethane foam cabinet insulation</li> <li>• Hydrocarbon, natural refrigerant (R290)</li> <li>• Pharmacy refrigerator/freezer toolkit and temperature logs</li> </ul>

**SPEC SHEET PH-ABT-HC-SSP-49**

<b>Product Name:</b>	PH-ABT-HC-SSP-49
<b>Features:</b>	<ul style="list-style-type: none"><li>• Two year parts and labor warranty, plus an additional three year compressor parts warranty</li><li>• Exterior dimensions: 54" W x 35" D x 84" H with casters; 77 ½" H without casters</li><li>• Cycle defrost</li><li>• 115V, 60 Hz, 4.5 Amps, 1/4 HP</li><li>• ETL listed</li><li>• Shipping weight: 450 lbs.</li></ul>