

SPEC SHEET PH-ABT-HC-SSP-49FA

Product Name:	PH-ABT-HC-SSP-49FA
Features:	<ul style="list-style-type: none"> • Premier Series Pharmacy Stainless Steel Freezer • 49 Cu. Ft. capacity • Two swing solid doors, self-closing, not reversible • Six ea. Adjustable Shelves • Digital display microprocessor temperature controller • Adjustable operating temperature: 5°F – -13°F (-15°C - -25°C) • Temperature Monitor Device - complies with the current CDC guidelines, 3 years certification of calibration, “buffered” probe in the product simulated solution, min/max memory, °F/°C switchable, field installable, and visual & audible temperature alarms • 2 probes; 1 in air and 1 in sample bottle • Probe attached to evaporator for defrost termination sensing • Sensor error alarm • Continuous product (simulated) temperature display • Remote alarm contacts • Forced draft air circulation • Keyed door lock • Magnetic door gaskets for positive seal • Integrated door handle • LED interior lights are safety shielded and switch controlled • Probe access port (3/4") • Stainless steel exterior • White-powder coated aluminum interior with 304 stainless steel floor • Easily cleaned smooth interior • Hydrocarbon natural refrigerant (R290)

SPEC SHEET PH-ABT-HC-SSP-49FA

Product Name:	PH-ABT-HC-SSP-49FA
Features:	<ul style="list-style-type: none"> • 2 year parts and labor warranty, plus an additional 3 year compressor parts warranty • Auto defrost • 4 swivel casters (5 3/4") • 115V/60 Hz • ETL listed • Exterior dimensions: 54" W x 35 3/4" D x 83 1/4" H with casters; 77 1/2" H without casters • 115V, 60 Hz, 9 Amps, 3/4 HP • ETL listed • Shipping weight: 495 lbs.