

# SPEC OF AT250A-9-N

STP-JX-Tip-AT250A-9-A1

GEB Robotic Tips for Agilent® / Agilent Bravo®.

Product details at-a-glance

GEB Cat#	Max Volume	Tip Color	Low Retention	Sterilized	Filtered (4.1*3.5)mm	Packaging Spec.
AT250A-9-N	250ul	Natural				96 tips/rack, 50 racks/case 4,800 tips/case 5 packs/case.
AT250A-9-NS	250ul	Natural		Sterilized		
AF250A-9-NS	250ul	Natural		Sterilized	Filtered	
AF250A-9-NS-LB	250ul	Natural	Low Retention	Sterilized	Filtered	
Case dimension L*W*H: 500*285*405 mm.			Gross Weight : 13.0 kg			

- ✧ 250µL Natural Tips
- ✧ 96 Tip/Hanging Tip Rack format
- ✧ Polypropylene
- ✧ Certified DNase/RNase free
- ✧ Nonpyrogenic
- ✧ Presterile/Nonsterile
- ✧ Filtered(Filtered With aerosol-resistant filters)/NonFiltered
- ✧ Endotoxin-Safe meets USP/FDA standards
- ✧ With USP 3127 (Heavy Metal Contamination)
- ✧ Compatibility: Agilent® / Agilent Bravo® workstations.

Drawing: (No tolerance Diameter tolerance:  $\pm 0.05\text{mm}$  , Height tolerance:  $\pm 0.25\text{mm}$ )

