

# SPEC OF AT300H-H2-X

STP-JX-Tip-AT300H-H2-A1

GEB Robotic Tips for Hamilton STAR®, STAR, STARlet, STARplus and Nimbus®. Laboratory Use Only.

Product details at-a-glance

GEB Cat#	Max Volume	Tip Color	Low Retention	Sterilized	Filtered (2.8*3.5)mm	Packaging Spec.
AUH000300	300ul	Black				96 tips/rack, 2,304 tips/case
ANH000300	300ul	Black	Low Retention		Filtered	
AUH001300	300ul	Black		Sterilized		
ANH001300	300ul	Black	Low Retention			
AUH101300	300ul	Black		Sterilized	Filtered	
ANH101300	300ul	Black	Low Retention	Sterilized	Filtered	
AT300H-H2-X	300ul	Black				2*96 tips/rack, 5 racks/pack, 10 packs/case 9600 tips/case. (Conductive black)
AT300H-H2-X-LB	300ul	Black	Low Retention			
AT300H-H2-XS	300ul	Black		Sterilized		
AT300H-H2-XS-LB	300ul	Black	Low Retention	Sterilized		
AF300H-H2-XS	300ul	Black		Sterilized	Filtered	
AF300H-H2-XS-LB	300ul	Black	Low Retention	Sterilized	Filtered	
Case dimension L*W*H: 405*285*235 mm.			Gross Weight : 11.6kg			

- ✧ 300µL Liquid Level Sensing Tips
- ✧ 96 Tip/Hanging Tip Rack format
- ✧ Double stack box packaging
- ✧ Polypropylene
- ✧ Certified DNase/RNase free
- ✧ Nonpyrogenic
- ✧ Presterile/Nonsterile
- ✧ Filtered(Filtered With aerosol-resistant filters)/NonFiltered
- ✧ Endotoxin-Safe meets USP/FDA standards
- ✧ Comply With USP 3127 (Heavy Metal Contamination) standards
- ✧ Compatibility: Hamilton STAR®, STAR, STARlet, STARplus and Nimbus® workstations.

Drawing: (No tolerance Diameter tolerance: ±0.03mm , Height tolerance: ±0.25mm)

