

# SPEC OF AT050H-H2-X

STP-JX-Tip-AT050H-H2-A1

GEB Robotic Tips for Hamilton STAR®, STAR, STARlet, STARplus and Nimbus®.

## Product details at-a-glance

GEB Cat#	Max Volume	Tip Color	Low Retention	Sterilized	Filtered (2.8*3.5)mm	Packaging Spec.
AUH000050	50ul	Black				96 tips/rack, 2,304 tips/case
ANH000050	50ul	Black	Low Retention		Filtered	
AUH001050	50ul	Black		Sterilized		
ANH001050	50ul	Black	Low Retention			
AUH101050	50ul	Black		Sterilized	Filtered	
ANH101050	50ul	Black	Low Retention	Sterilized	Filtered	
AT050H-H2-X	50ul	Black				2*96 tips/rack, 5 racks/pack, 10 packs/case 9600 tips/case. (Conductive black)
AT050H-H2-X-LB	50ul	Black	Low Retention			
AT050H-H2-XS	50ul	Black		Sterilized		
AT050H-H2-XS-LB	50ul	Black	Low Retention	Sterilized		
AF050H-H2-XS	50ul	Black		Sterilized	Filtered	
AF050H-H2-XS-LB	50ul	Black	Low Retention	Sterilized	Filtered	
Case dimension L*W*H: 465*295*455 mm.			Gross Weight : 11.6 kg			

- ✧ 50µL Liquid Level Sensing Tips
- ✧ 96 Tip/Hanging Tip Rack format
- ✧ Double stack box packaging
- ✧ Polypropylene
- ✧ Certified DNase/RNase free
- ✧ Nonpyrogenic
- ✧ Presterile/Nonsterile
- ✧ Filtered(Filtered With aerosol-resistant filters)/NonFiltered
- ✧ Endotoxin-Safe meets USP/FDA standards
- ✧ Comply With USP 3127 (Heavy Metal Contamination) standards
- ✧ Compatibility: Hamilton STAR®, STAR, STARlet, STARplus and Nimbus®workstations.

Drawing: (No tolerance Diameter tolerance: ±0.03mm , Height tolerance: ±0.25mm)

