

HandyStep®

Repeating Pipettes



BRAND® HandyStep® touch/touch S

Choose the HandyStep® touch for your repetitive pipetting tasks—dispensing, auto-dispensing and the ability to store favorites. The HandyStep® touch S adds such advanced features like multi-aspiration, sequential dispensing and titration. Paired with the new PD-Tip™ // Dispenser Tips with tip size coding, they make the perfect team for easy pipetting.

- **Works with All Standard Tips:** The HandyStep® touch is compatible with BRAND® PD-Tip™ // Precision Dispenser Tips, Eppendorf® Combitips advanced®, and more.
- **Recognizes Tips Automatically:** The pipette automatically recognizes the size of encoded BRAND® PD-Tip™ // Precision Dispenser Tips upon tip insertion and displays the tip's volume, default dispensing volume, and default number of dispensing steps.
- **Reduces Repetitive Motion by up to 97%:** The instrument features a unique Auto-Dispense operating mode in which the pipette "learns" the operator's dispensing rhythm and dispenses at that rhythm as long as the control button is depressed – no programming needed!

Use the tips you want

The only touch screen repeating pipette compatible with all standard syringe tips, the HandyStep® touch accepts:

- BRAND® PD-Tip™ // Precision Dispenser Tips
- Rainin® Encode™ syringe tips
- Eppendorf® Combitips advanced®
- Fisherbrand® Dispenser Tips
- VWR® Combi-Syringes
- Gilson® Repet-Tips

Description	Cat. No.	2022 List Price
HandyStep® S	705110	\$426.00
HandyStep® touch	705200	1,109.00
HandyStep® touch S	705210	1,374.00

BRAND® HandyStep® S

The HandyStep® S repeating pipette provides unprecedented versatility in a purely mechanical repeating pipette. In conjunction with BRAND® PD-Tip™ // Precision Dispenser Tips, the HandyStep® S offers nearly sixty different volume settings with accuracy, precision and dependability.

- **Comfortable:** The slim design of the HandyStep® S is lightweight, well-balanced, and provides a convenient control layout for smooth dispensing and one-handed volume selection.
- **Reliable:** Robust construction, inspired design and innovative polymers deliver a dependable instrument.
- **Ready-To-Go:** A purely mechanical design with no batteries to charge or replace means the HandyStep® S is ready when you are. A convenient integrated volume table provides a quick reference for dispensed volumes.
- **Flexible:** Ten sizes of BRAND® PD-Tip™ // Precision Dispenser Tips allow the HandyStep® S to dispense 59 different volumes. Open design allows use of third party tips—no more being locked into a sole source!
- **Contact-Free:** A tip ejection system allows the user to remove tips without handling.

Description	Cat. No.	2022 List Price
HandyStep® S	705110	\$426.00
HandyStep® touch	705200	1,109.00
HandyStep® touch S	705210	1,374.00



PD-Tip™ // Syringe Tips

For Repeating Pipettes

The Perfect Tips for Your Repeating Pipette.

BRAND® encoded PD-Tip™ // syringe tips are manufactured from high-quality virgin plastics for use with all standard repeating pipettes in a wide variety of applications.

- **Compatible with Most Pipettes:** BRAND® PD-Tip™ // Precision Dispenser Tips are compatible with repeating pipettes using standard tips. BRAND® PD-Tip™ // Precision Dispenser Tips are encoded for automatic recognition by the BRAND® HandyStep® touch repeating pipette, the Rainin® AutoRep™ E pipette, and the Gilson® Repetman®.
- **Reliable:** A precision seal between the piston and cylinder ensure smooth operation, increased accuracy, and increased precision. Extended piston design minimizes dead volume and maximizes fluid utilization.
- **Manufactured from High Quality Materials:** PD-Tip™ // syringe tips are made entirely from high-quality virgin polyethylene and polypropylene. (Piston of 0.1 mL tip is liquid crystal polymer-LCP).
- **BIO-CERT® PD-Tip™ // Precision Dispenser Tips** are sterile, individually wrapped, and certified free of DNA, RNase, endotoxins, and ATP to meet the most stringent sample integrity requirements.
- **Batch Certificate of Performance** which documents tip accuracy and precision testing is supplied with each pack of BRAND® PD-Tip™ // Precision Dispenser Tips.



PD-Tip™ // Precision Dispenser Tips BIO-CERT® sterile, individually wrapped

Capacity	Pack of	Item #	2022 List Price
0.1 mL	100	705730	\$182.80
0.5 mL	100	705732	\$182.80
1 mL	100	705734	\$182.80
1.25 mL	100	705736	\$182.80
2.5 mL	100	705738	\$182.80
5 mL	100	705740	\$182.80
10 mL	100	705742	\$182.80
12.5 mL	100	705744	\$182.80
25 mL (one adapter included)	25	705746	\$56.00
50 mL (one adapter included)	25	705748	\$56.00
Adapter for 25 mL and 50 mL , BIO-CERT® sterile, PP, autoclavable			
25 mL and 50 mL Adapter	5	702399	\$83.40

*The Eppendorf® Repeater® Plus, Repeater® M4, Repeater® Pro, Stream, and XStream are designed to preclude their operation with any tips other than Eppendorf® Combitips® advanced.

PD-Tip™ // Precision Dispenser Tips, non-sterile

Capacity	Pack of	Item #	2022 List Price
0.1 mL	100	705700	\$133.40
0.5 mL	100	705702	\$133.40
1 mL	100	705704	\$133.40
1.25 mL	100	705706	\$133.40
2.5 mL	100	705708	\$133.40
5 mL	100	705710	\$133.40
10 mL	100	705712	\$133.40
12.5 mL	100	705714	\$133.40
25 mL (one adapter included)	50	705716	\$90.40
50 mL (one adapter included)	25	705718	\$45.20
Multi-pack (20 each 0.5, 1, 1.25, 2.5, 5, 10 and 12.5 mL)			
	140	705720	\$169.60
Adapter for 25 mL and 50 mL , non-sterile, PP, autoclavable			
25 mL and 50 mL	10	702398	\$99.60



Lab Rats Trust BrandTech!

HandyStep_BT_Fly-2022

BRANDTECH®
SCIENTIFIC, INC.

MIDSCI™
Laboratory Equipment & Supplies

ORDER TODAY! midsci.com • custserv@midsci.com • 800.227.9997 • f: 636.225.9998