

SPEC OF AT800R-B/9

STP-JX-Tip-AT800R-B/9-A1

GEB Robotic Tips for Rainin® can be used in proteomics, drug development, genomics and other fields.

Product details at-a-glance

GEB Cat#	Max Volume	Tip Color	Low Retention	Sterilized	Filtered (4.5*3.5)mm	Packaging Spec.
AT800R-9-X	800ul	Black				96 tips/rack, 10 racks/pack, 5 packs/case. 50 racks/case, 4,800 tips/case.
AT800R-9-XS	800ul	Black		Sterilized		
AF800R-9-XS	800ul	Black		Sterilized	Filtered	
AF800R-9-XS-LB	800ul	Black	Low Retention	Sterilized	Filtered	
Case dimension L*W*H::480* 345*280mm			Gross Weight : 10.0 kg			

- ✧ 800µL Liquid Level Sensing Tips
- ✧ 96 Tip/Hanging Tip Rack format
- ✧ Polypropylene
- ✧ Certified DNase/RNase free
- ✧ Nonpyrogenic
- ✧ Presterile/Nonsterile
- ✧ Filtered(Filtered With aerosol-resistant filters)/NonFiltered
- ✧ Endotoxin-Safe meets USP/FDA standards
- ✧ Comply With USP 3127 (Heavy Metal Contamination) standards
- ✧ Compatibility: Rainin® liquid workstations.

Drawing: (No tolerance Diameter tolerance : ±0.05mm , Height tolerance : ±0.25mm)

