



- Digital speed control up to 300 rpm
- Extremely compact, perfect for low to medium capacity applications
- Instantly exchange flask clamps with MAGic Clamp™ accessories
- Horizontally circular 19 mm orbit for aeration & mixing

The PR1MA™ Orbi-Shaker™ Small is designed for smaller capacity applications and laboratories where bench space is at a premium. The shaker includes a 9.5 x 8 in. platform that mixes at a 19 mm orbit, perfect for thorough mixing and aeration in Erlenmeyer flasks.

An optional MAGic Clamp™ platform is available for use with Erlenmeyer flasks and test tube racks. The unique magnetic attachment method is the easiest way to instantly change between flask clamps of different sizes.

SPECIFICATIONS

SPEED	Variable, 50 to 300 rpm
MAX CAPACITY	4 x 1 L (9 lbs / 4 kg)
DIMENSIONS	13.1 x 11.74 x 4.7 in. 33 x 30 x 12 cm
WEIGHT	12 lbs (5.5 kg)
ELECTRICAL	Universal Voltage, 100-240 V, 50-60 Hz 250 W
WARRANTY	2-Year Warranty



H1000-MR-CMB

ORDERING INFORMATION

ITEM #	DESCRIPTION	QUANTITY
OS-300	Variable Speed, 9.5 x 8 in. w/ Rubber Mat Platform	1/ea
H1000-MR-CMB	MAGic Clamp™ Combination Clamp	1/ea
IS-PLAT	MAGic Clamp™ Platform for OS-300	1/ea