



- Powerful motor for INSTANT vortexing
- Continuous or “touch” operation
- Q-Drive™ dynamic balancing system
- Variable speed from 200 to 3,200 rpm
- Wide variety of accessories

The PR1MA™ Vortex Mixer sets the new standard in vortex mixing. With advanced Q-Drive™ technology, the PR1MA™ Vortex Mixer provides smooth instant vortexing of even the largest sample sizes, including full 50 mL tubes. The unique counter-balance system creates maximum vortexing action while minimizing noise and excessive vibration. This results in quieter, more efficient, and longer-lasting operation.

Other features include “no-walk” suction feet, cold room/incubator compatibility, and spill-proof electronic controls.



Unique multi-head allows for standard vortexing PLUS holds up to 8 microtubes

SPECIFICATIONS

SPEED	200 to 3200 rpm
ORBIT	3 mm
DIMENSIONS	5 x 6.3 x 6.7 in. 13 x 16 x 17 cm
WEIGHT	8.4 lbs 3.8 kg
ELECTRICAL	115 V, 60 Hz or 230 V, 50 Hz, 180 W
WARRANTY	2-Year Warranty

ORDERING INFORMATION

ITEM #	DESCRIPTION	QUANTITY
VM-3200	Adjustable Vortex Mixer, 3 mm Orbit, Max Speed 3,200 rpm	1/ea
VM-COMBO	Combo Head Attachment, 1.5 mL, 0.5 mL Microtubes & 1 Microplate	1/ea
VM-FLAT	Flat Head Attachment, 3 in.	1/ea
VM-H15	Horizontal Head, 12 x 1.5 mL Microtubes	1/ea
VM-H150	Horizontal Head, 4 x 15 mL Microtubes	1/ea
VM-H500	Horizontal Head, 2 x 50 mL Microtubes	1/ea