




TX-1000, adapters, buckets, and bio lids for the smaller units to accommodate 400-450 mL Bioreactor Tubes.



TX-1000 4 x 1000mL



75007302



75007309



75003001

Ausschwingrotor / Swing-Out Rotor		TX-1000 Rotor Cross	75003017
Becher / Bucket 75003001	Kappen / Bucket Lid 75007309	Adapter 75007302	Kapazität / Capacity 4 x 87450

*This setup will work in any of the centrifuges listed below.*

Centrifuge	
Sorvall X4R/X4FR Pro	ST 40 SL 40R
Sorvall X4/X4F Pro	SL 40 SL40FR
Sorvall ST4R/ST4FR Plus	SL 40F
Sorvall ST4/ST4F Plus	
SL4R/SL4FR Plus	Multifuge X4R Pro/X4FR Pro
SL4/SL4F Plus	Multifuge X4 Pro/X4F Pro
	Megafuge ST4R/ST4FR Plus
	Megafuge ST4/ST4F Plus
Legend XTR	Multifuge X3R
Legend XT	Multifuge X3
Legend XFR	Multifuge X3FR
Legend XF ST 40R	Multifuge X3F
	Multifuge 40R
	Megafuge 40



For the 600 mL bottle, it would be these two rotors along with the buckets and adapters:



HAEMAFlex 6 Swinging Bucket Rotor



HAEMAFlex 8 Swinging Bucket Rotor

Ausschwingrotor / Swing-Out Rotor		HAEMAFlex 6	75003861
		HAEMAFlex 8	75003881
Becher / Bucket 75004871	Kappen / Bucket Lid ---	Adapter 75004872	Kapazität / Capacity 6 / 8 x 87600

Centrifuge
Sorvall BP 8
Sorvall BP 8 Heavy Duty
Sorvall BP 8 Heavy Duty Water Cooled
Cryofuge 8
Cryofuge 8 Heavy Duty
Cryofuge 8 Heavy Duty Water Cooled
Sorvall BIOS 8
Sorvall BIOS 8 Heavy Duty

*\*The bucket (75004871) for the Larger capacity centrifuges will be available in North America in 2024.*

*Please contact us with questions or to inquire about availability.*