


Round Bottom Snap Cap Tubes With Attached Caps

- Microtube convenience in larger capacities
- Attached caps seal at the rim with a "click"
- Sterile, easy-open bags of 100 tubes
- Molded graduations and frosted marking area



MATERIAL:	Polypropylene
TEMPERATURE RANGE:	-86°C - 121°C
CENTRIFUGATION RATING:	Up to 4,000×g
DIMENSIONS:	5mL - 12×55mm 7mL - 13×62mm 10mL - 14×82mm

Stop fumbling with caps for your tubes!

MTC Bio's new snap cap culture/centrifuge tubes adopt the conveniences of popular 1.5mL microcentrifuge tubes in higher capacities of 5mL, 7mL, and 10mL. Forget the cross-contamination concerns of mismatched caps!

These durable polypropylene tubes feature molded graduations on all sizes. The smaller 5mL and 7mL tubes also have a frosted marking area. Tubes are supplied in sterile, easy-open bags of 100 and are autoclavable.



CAT. NO.	DESCRIPTION	QTY
C0905	5mL Sterile, w/ Attached Snap Cap, Molded Graduations, Frosted Marking Area	1000 / case
C0907	7mL Sterile, w/ Attached Snap Cap, Molded Graduations, Frosted Marking Area	1000 / case
C0910	10mL Sterile, w/ Attached Snap Cap, Molded Graduations	1000 / case