

## Haier Biomedical Equipment Specifications

Product	ULT freezer model	Rack model	Rack and box description	Freezer Rack Quantity	Freezer Box Quantity	Cryo Tube Quantity (2 mL)
<b>Racks for ULT Freezers</b>	DW-86L338J /DW-86L338JA	DCJ-54-A (SUS304)	2 inch stainless steel side-access freezer rack 5*4 (2 inch freezer box*20 units)	6	120	12,000
		DCJ-54-B (SUS201/304) (Middle, 2 layers)	2 inch stainless steel sliding-drawer freezer rack 5*4 (2 inch freezer box*20 units)	6	120	12,000
		DCJ-44-A (SUS304) (Top and bottom layers)	2 inch stainless steel side-access freezer rack 4*4 (2 inch freezer box*16 units)	6	96	9,600
		DCJ-44-B (SUS201/304) (Top and bottom layers)	2 inch stainless steel sliding-drawer freezer rack 4*4 (2 inch freezer box*16 units)	6	96	9,600
	DW-86L388J/418S /418ST/419	DCJ-55-A (SUS304)	2 inch stainless steel side-access freezer rack 5*5 (2 inch freezer box*25 units)	12	300	30,000
		DCJ-55-B (SUS201/304)	2 inch stainless steel sliding-drawer freezer rack 5*5 (2 inch freezer box*25 units)	12	300	30,000
	DW-86L486E	DCJ-54-A (SUS304)	2 inch stainless steel side-access freezer rack 5*4 (2 inch freezer box*20 units)	16	320	32,000
		DCJ-54-B (SUS201/304)	2 inch stainless steel sliding-drawer freezer rack 5*4 (2 inch freezer box*20 units)	16	320	32,000
	DW-86L490J/490JA	DCJ-54-A (SUS304) (Top 2 layers)	2 inch stainless steel side-access freezer rack 5*4 (2 inch freezer box*20 units)	8	160	16,000
		DCJ-54-B (SUS201/304) (Top 2 layers)	2 inch stainless steel sliding-drawer freezer rack 5*4 (2 inch freezer box*20 units)	8	160	16,000
		DCJ-44-A (SUS304) (Bottom 2 layers)	2 inch stainless steel side-access freezer rack 4*4 (2 inch freezer box*16 units)	8	128	12,800
		DCJ-44-B (SUS201/304) (Bottom 2 layers)	2 inch stainless steel sliding-drawer freezer rack 4*4 (2 inch freezer box*16 units)	8	128	12,800
	DW-86L628E	DCJ-54-A (SUS304)	2 inch stainless steel side-access freezer rack 5*4 (2 inch freezer box*20 units)	20	400	40,000
		DCJ-54-B (SUS201/304)	2 inch stainless steel sliding-drawer freezer rack 5*4 (2 inch freezer box*20 units)	20	400	40,000

Product	ULT freezer model	Rack model	Rack and box description	Freezer Rack Quantity	Freezer Box Quantity	Cryo Tube Quantity (2 mL)
<b>Racks for ULT Freezers</b>	DW-86L578J/578JA /578S/578ST/579 /579BP/579BPT	DCJ-55-A (SUS304)	2 inch stainless steel side-access freezer rack 5*5 (2 inch freezer box*25 units)	16	400	40,000
		DCJ-55-B (SUS201/304)	2 inch stainless steel sliding-drawer freezer rack 5*5 (2 inch freezer box*25 units)	16	400	40,000
	DW-86L728J/728JA /728S/728ST/728SAT /729/729BP/729BPT	DCJ-55-A (SUS304)	2 inch stainless steel side-access freezer rack 5*5 (2 inch freezer box*25 units)	20	500	50,000
		DCJ-55-B (SUS201/304)	2 inch stainless steel sliding-drawer freezer rack 5*5 (2 inch freezer box*25 units)	20	500	50,000
	DW-86L828J/828JA /828S/828ST/829 /829BP/829BPT/828W	DCJ-55-A (SUS304)	2 inch stainless steel side-access freezer rack 5*5 (2 inch freezer box*25 units)	24	600	60,000
		DCJ-55-B (SUS201/304)	2 inch stainless steel sliding-drawer freezer rack 5*5 (2 inch freezer box*25 units)	24	600	60,000
	DW-86L959BP /959BPT/959W	DCJ-55-A (SUS304)	2 inch stainless steel side-access freezer rack 5*5 (2 inch freezer box*25 units)	28	700	70,000
		DCJ-55-B (SUS201/304)	2 inch stainless steel sliding-drawer freezer rack 5*5 (2 inch freezer box*25 units)	28	700	70,000
	DW-86L100J	DCJ-53-B(SUS201/304)	2 inch stainless steel sliding-drawer freezer rack 5*3 (2 inch freezer box*15 units)	4	60	6,000
	DW-86W420J/420JA	DCJ-10-A (SUS304)	2 inch chest-type freezer rack (2 inch freezer box*10 units)	27	270	27,000
	DW-86W100J	DCJ-08-A (SUS304)	2 inch chest-type freezer rack (2 inch freezer box*8 units)	9	72	7,200
	DW-150W200	DCJ-10-A (SUS304)	2 inch chest-type freezer rack (2 inch freezer box*10 units)	12	120	12,000
	DW-150W209	DCJ-10-A (SUS304)	2 inch chest-type freezer rack (2 inch freezer box*10 units)	12	120	12,000
		DCJ-12-A (SUS304)	2 inch chest-type freezer rack (2 inch freezer box*12 units)	12	144	14,400
	DW-86L51J	DW-86L51J	2 inch stainless steel sliding-drawer freezer rack 5*3 (2 inch freezer box*15 units)	2	30	3,000