



Laboratory Equipment & Supplies

Offering the best values on ELISA instruments

*Free 8-channel Pipette w/purchase



Solutions for ELISA reading and washing



- Intuitive, 7 inch color touch screen
- Filter based system, wavelengths from 340 to 750nm
- Absorbance range: 0.0 to 4.000 Abs
- 8 channel vertical optical path with zero dispersion

- Precision control of aspiration needle position
 - Accurate to 0.1mm
 - Minimal residual liquid, $\leq 1\mu\text{l/well}$
- Auto rinse function eliminates cross contamination
- Up to 100 wash protocols can be customized and stored
- Liquid level sensors in wash and waste bottles

MR9600-MID	SmartReader™96 Microplate Reader	\$5,395.00
MR9600-T	SmartReader™96 With Temperature control	\$5,750.00


MW9600	SmartWasher™96 Microplate Washer	\$4,295.00
--------	---	-------------------

* Choose 8 and 12 channel manifold, sold separately



Free 8-channel pipette with purchase of any of these instruments. Use promo code **ELISA8** (Value: \$389.00)

Solutions for ELISA plate mixing and incubating



- Two models: Heat/Cool or Heat only
- Temperature control from 0° to 100°C
- Exchangeable blocks, from 0.2ml to 50ml

- Two models: two-plate or four-plate format
- Heated lid prevents condensation
- Temperature control up to 70°C
- Accepts a variety of microplates

MTH-1500	MultiTherm™ Shaker with heating only	\$1,999.00
MTHC-1500	MultiTherm™ Shaker with heating and cooling	\$2,629.00
MTB-MP	Block for on Micro Plate	\$385.00

*Additional blocks available

IM6002	Incu-Mixer MP2™ Heated Plate Vortexer	\$1,635.00
IM6004	Incu-Mixer MP4™ Heated Plate Vortexer	\$2,249.00

*Special offers valid through December 31, 2020

ORDER TODAY! midsci.com • custserv@midsci.com • 800.227.9997 • f: 636.225.9998