




- Twice the capacity of traditional mini centrifuges
- Near silent operation
- Starts and stops with closing/opening of the lid

Spin microtubes and PCR strips simultaneously without changing rotors!

Unlike traditional mini centrifuges, the PR1MA™ 12 eliminates the need to change rotors when switching between microtubes and PCR strips. The included, unique COMBI-Rotor™ is all that is required for running 12 microtubes and 4 PCR strips simultaneously.



Compact, low profile design



Unique COMBI-Rotor™

SPECIFICATIONS

| | |
|---------------------|---|
| SPEED | 5,500 rpm 2,000 x g |
| MAX CAPACITY | 12 x 1.5/2.0 mL Microtubes and 32 x 0.2 mL PCR Tubes (or 4 PCR Strip Tubes) |
| DIMENSIONS | 5.5 x 7.9 x 4.4 in. |
| WEIGHT | 2.1 lbs / 5 kg |
| WARRANTY | 2-Year Warranty |

ORDERING INFORMATION

| ITEM # | DESCRIPTION | QUANTITY |
|--------|--|----------|
| MF12 | 12 Mini, Capacity 12 x 1.5 mL Tubes, Max Speed 5,500 rpm | 1/ea |