

PR1MA™ 8 Mini Microcentrifuge

PR1MA™



- Clik-n-Lok™ Instant rotor attachment
- Eight position microtube rotor
- Decelerates in just 1 second
- Storage compartment for PCR rotor and 0.5 mL adapters
- Four position PCR strip rotor

The PR1MA™ Mini centrifuge is a personal centrifuge that virtually fits in the palm of your hand. Ideal for quick spin downs of microtubes and PCR tubes, the PR1MA Mini centrifuge is extremely easy to use. Simply close the lid and the rotor instantly reaches 6,000 rpm. Open the lid, and the improved braking system quickly and smoothly decelerates the rotor to a complete stop in just one second.

An eight-position microtube rotor is included and comes installed. Also included is a second rotor with increased capacity for four 0.2 mL PCR strips (or 32 x 0.2 mL tubes) which is stored in a convenient compartment on the bottom of the centrifuge. This compartment is also used for securing the included tube adapters for 0.5 mL tubes.

A new Clik-n-Lok rotor attachment system eliminates the requirement for set screws and tools during rotor exchange. Simply remove and attach rotors by hand, pulling or pressing the rotor onto the motor shaft.



Includes Both Rotors



Storage Compartment

SPECIFICATIONS

SPEED	6000 rpm / 2000 x g
CAPACITY	8 x 1.5/2.0 mL Tubes 32 x 0.2 mL Tubes 4 x PCR Strips
DIMENSIONS	5.7 x 6.9 x 4.5 in. 14.5 x 17.8 x 11.4 cm
WEIGHT	2.1 lbs / 5 kg
ELECTRICAL	100 to 240V, 50-60 Hz

ORDERING INFORMATION

ITEM #	DESCRIPTION	QUANTITY
MF8-B	PR1MA™ Mini Microcentrifuge, Blue Lid	1/ea
MF8-C	PR1MA™ Mini Microcentrifuge, Clear Lid	1/ea
MF8-G	PR1MA™ Mini Microcentrifuge, Green Lid	1/ea
MF8-P	PR1MA™ Mini Microcentrifuge, Purple Lid	1/ea
MF8-R	PR1MA™ Mini Microcentrifuge, Red Lid	1/ea

All Include microtube rotor, PCR rotor, 0.5 mL & 0.2 mL adapters

MIDSCI™
Laboratory Equipment & Supplies

ORDER TODAY! midsci.com • custserv@midsci.com • 800.227.9997 • f: 636.225.9998