

PR1MA™ Mini Incubated Shaker

PR1MA™



- Redesigned to accept two 1 L flasks
- Increased temperature control up to 70°C
- Instant clamp exchange with MAGIC Clamp™ system

The PR1MA™ Mini Incubated Shaker is the most compact shaking incubator in its class. Despite the remarkably small footprint, a large 11.5 x 9.5 in. work space accepts several platforms for use with a variety of vessels.

The standard platform (included) features a non-slip rubber-coated surface ideal for TC flasks, petri dishes, and staining trays.

The convenient MAGIC Clamp™ Platform is also available for use with Erlenmeyer flasks and test tube racks. The unique, magnetic attachment method is the easiest way to instantly change between flask clamps of different sizes.



SPECIFICATIONS

| | |
|----------------------------|---|
| SHAKING SPEED RANGE | 30-300 rpm |
| TEMP. RANGE | Ambient + 5°C to 70°C (0.1°C increments) |
| DIMENSIONS | 11 x 15.7 x 10.5 in. (height 25 in. with lid open) |
| WEIGHT | 25 lbs |
| ELECTRICAL | 100 to 240V, 50-60 Hz |
| WARRANTY | 2-Year Warranty |



ORDERING INFORMATION

| ITEM # | DESCRIPTION | QUANTITY |
|--------------|--|----------|
| IS-300 | Mini Incubated Shaker, Ambient + 5°C to 70°C, Max Speed 300 rpm, 9.5 x 11.5 in. Platform | 1/ea |
| IS-PLAT | MAGIC Clamp™ Platform for PR1MA™ Mini Incubated Shaker | 1/ea |
| H1000-MR-CMB | MAGIC Clamp™ Combination Clamp | 1/ea |

MIDSCI™
Laboratory Equipment & Supplies

ORDER TODAY! midsci.com • custserv@midsci.com • 800.227.9997 • f: 636.225.9998