

PR1MA™ Mini Vortexer

PR1MA™



- Powerful vortexing for tubes up to 50 mL
- Pressure activated operation
- "Mini-sized", less than 4 inches wide
- Four position PCR strip rotor

Powerful is the best word to describe the PR1MA™ Mini Vortexer. With a 4 mm orbit and fixed speed of 2,800 rpm, it instantly vortexes even the largest samples, including nearly full 50 mL tubes. Despite its powerful motor, the Mini Vortexer has a footprint of less than 4 x 4 inches, allowing it to fit on even the most crowded lab bench.

Pressure-activated operation means that the motor is automatically activated with slight pressure on the center cup head. Operation is smooth, quiet, and efficient. When pressure is removed operation immediately ceases. The unique design prevents liquids from entering the housing, prolonging motor life.



Extremely Compact



Insta-Touch™ Pressure Activation

SPECIFICATIONS

SPEED	2,800 rpm
ORBIT	4 mm
DIMENSIONS	3.7 x 3.9 x 2.6 in.
WEIGHT	1 lb.
ELECTRICAL	Universal Voltage, 100-240 V, 50-60 Hz, 0.7 A
WARRANTY	2-Year Warranty

ORDERING INFORMATION

ITEM #	DESCRIPTION	QUANTITY
VM1-B	Mini Vortexer, Blue, 4 mm Orbit, Max Speed 2,800 rpm	1/ea
VM1-G	Mini Vortexer, Green, 4 mm Orbit, Max Speed 2,800	1/ea
VM1-GR	Mini Vortexer, Gray, 4 mm Orbit, Max Speed 2,800	1/ea
VM1-P	Mini Vortexer, Purple, 4 mm Orbit, Max Speed 2,800	1/ea
VM1-R	Mini Vortexer, Red, 4 mm Orbit, Max Speed 2,800	1/ea

MIDSCI™
Laboratory Equipment & Supplies

ORDER TODAY! midsci.com • custserv@midsci.com • 800.227.9997 • f: 636.225.9998