

COMFORT AND DEXTERITY FOR EVERYDAY HAND PROTECTION.



TF-95RB

**POWDER-FREE ECOTEK BIODEGRADABLE
NITRILE EXAMINATION GLOVES***


VALUE AND PERFORMANCE HAND PROTECTION FOR LAB, LIFE SCIENCE AND INDUSTRIAL ENVIRONMENTS

- EcoTek® Biodegradability Technology accelerates landfill breakdown without any performance loss*
- Silicone-free for critical environments
- Exceeds 1.5 AQL laboratory standards

The right weight
for resiliency and
dexterity

Fully textured
surface for wet
and dry grip

Tear-resistant
beaded cuff for
easy donning



TF **TF-95RB**

WEIGHT	LENGTH	THICKNESS	FEATURES	COLOR
5.0g (large size)	9.5"	3.1 mil	EcoTek®	ROYAL BLUE

ECOTEK





TF-95RB

POWDER-FREE ECOTEK BIODEGRADABLE
NITRILE EXAMINATION GLOVES*

PRODUCT DESCRIPTION

SPECIFICATIONS

PRODUCT	TrueForm®
GRADE	Powder-Free Exam
MATERIAL	Nitrile
THICKNESS FINGER	4.7mil 0.12mm
PALM	3.1mi 0.08mm
COLOR	Royal Blue
SURFACE	Textured
LENGTH	9.5" 242mm
WEIGHT	5.0g
FIT TYPE	Ambidextrous
CUFF TYPE	Beaded
WATERTIGHT AQL	1.5
ORIGIN COUNTRY	Malaysia

PHYSICAL PROPERTIES

	PRE-AGING	POST-AGING
TENSILE STRENGTH:	14 MPA	14 MPA
ELONGATION:	500%	400%

INTEGRATED TECHNOLOGIES


ECOTEK Biodegradability Technology: 92.6% biodegradable in 2.5 years. ASTM D5526-12 test results do not predict future biodegradation.*



* Specifications based on size Large glove

STANDARDS & CERTIFICATIONS

EN ISO 21420:2020	EU 2016/425
EN 455-2000-2015	ANSI/ISEA 105-2016
ASTM D6319-19	CAT III Complex Design
ASTM D6978-05 (2019)	FDA 21 CFR 177.2600
ASTM F1671-13	
ASTM F739-12	
Reach EU 1907/2006	



Please visit swsafety.com/symbols for information on standards used on this document.

ORDERING

ITEM NO.	SIZE	PALM WIDTH	PACKAGING
N301341	XS 5½ – 6	3.0" 75mm	100 gloves per box 10 boxes per case
N301342	S 6½ – 7	3.3" 85mm	100 gloves per box 10 boxes per case
N301343	M 7½ – 8	3.7" 95mm	100 gloves per box 10 boxes per case
N301344	L 8½ – 9	4.1" 105mm	100 gloves per box 10 boxes per case
N301345	XL 9½ – 10	4.5" 115mm	100 gloves per box 10 boxes per case
N301346	2XL 10½ – 11	4.9" 125mm	100 gloves per box 10 boxes per case

APPLICATIONS

- INDUSTRIAL**
 - Food Handling and Packaging
 - Janitorial Maintenance
 - Painting
- MEDICAL**
 - Chemical Handling
 - Dental
 - Veterinary
 - Tattooing and Cosmetics

Keep out of sunlight. Store in a cool, dry place.
Keep away from sources of ozone or ignition.



R_N30134X1121

