



ENERGY STAR CERTIFIED

Lab Grade Refrigerators and Freezers

Haier - DW-86L728BPST : DW-86L728BPST

Specifications

ENERGY STAR Unique ID:	2470861
ENERGY STAR Partner:	Qingdao Haier Biomedical Co.,Ltd
Brand Name:	Haier
Model Name:	DW-86L728BPST
Model Number:	DW-86L728BPST
Product Type:	Ultra-Low Temperature Freezer
Total Volume (cubic feet):	25.42
Ultra-Low Temperature Freezer Energy Consumption (kWh/day/cu-ft):	0.35
Peak Temperature Variance (°C):	36.3
Product Form Factor:	Upright
Total Number of Doors:	1
Number of glass doors:	0
Number of solid doors:	1
Door Opening Orientation:	Swinging, Left
Number of Drawers or Shelves:	3
Number of Internal Ports and/or Access Holes:	4
Depth (in):	39.29
Height:	77.95
Width:	41.18
Defrost Type:	Manual
Refrigerant with GWP:	R-1150 (GWP:4), R-600a (GWP:3)
Installed Accessories:	N/A
Date Available on Market:	1682222400000
Date Certified:	1682568000000
Markets:	United States, Canada
ENERGY STAR Certified:	Yes

