



ENERGY STAR CERTIFIED

# Lab Grade Refrigerators and Freezers

## Haier - DW-86L578BPST : DW-86L578BPST

### Specifications

|  |                                  |
|--|----------------------------------|
| <b>ENERGY STAR Unique ID:</b>  | 2470720                          |
| <b>ENERGY STAR Partner:</b>  | Qingdao Haier Biomedical Co.,Ltd |
| <b>Brand Name:</b>   | Haier                            |
| <b>Model Name:</b>   | DW-86L578BPST                    |
| <b>Model Number:</b>   | DW-86L578BPST                    |
| <b>Product Type:</b>   | Ultra-Low Temperature Freezer    |
| <b>Total Volume (cubic feet):</b>  | 20.52                            |
| <b>Ultra-Low Temperature Freezer Energy Consumption (kWh/day/cu-ft):</b> | 0.41                             |
| <b>Peak Temperature Variance (°C):</b>                                   | 44.4                             |
| <b>Product Form Factor:</b>  | Upright                          |
| <b>Total Number of Doors:</b>  | 1                                |
| <b>Number of glass doors:</b>  | 0                                |
| <b>Number of solid doors:</b>  | 1                                |
| <b>Door Opening Orientation:</b>   | Swinging, Left                   |
| <b>Number of Drawers or Shelves:</b>                                     | 3                                |
| <b>Number of Internal Ports and/or Access Holes:</b>                     | 4                                |
| <b>Depth (in):</b>   | 39.29                            |
| <b>Height:</b>   | 77.95                            |
| <b>Width:</b>  | 35.24                            |
| <b>Defrost Type:</b>   | Manual                           |
| <b>Refrigerant with GWP:</b>   | R-1150 (GWP:4), R-600a (GWP:3)   |
| <b>Installed Accessories:</b>  | N/A                              |
| <b>Date Available on Market:</b>   | 1682222400000                    |
| <b>Date Certified:</b>   | 1682308800000                    |
| <b>Markets:</b>  | United States, Canada            |
| <b>ENERGY STAR Certified:</b>  | Yes                              |

