



ThermalSeal®

An easier to remove film

Application:

PCR and general storage


Temperature Range:

-40°C to +120°C

CATALOG #	DESCRIPTION	QUANTITY
100-THER-PLT	ThermalSeal	100
STR-THER-PLT	ThermalSeal, Sterile	100

ThermalSeal features a moderate strength adhesive for ease of removal.

- Certified DNase- & RNase- free
- Only recommended for structurally rigid polypropylene PCR plates that don't warp
- A less rigid polypropylene film that minimizes aerosols upon removal
- Features good chemical resistance



"For a stronger more secure seal, check out eXTReme Seal, Large"

