

AlumaSeal® II™

Easy-to-pierce for direct sample recovery

Application:

PCR, long-term storage, light sensitive assays and robotics


Temperature Range:

-80°C to +120°C

CATALOG #	DESCRIPTION	QUANTITY
AF-100	AlumaSeal II	100
AFS-25	AlumaSeal II, Sterile	50

AlumaSeal II is a soft non-permeable aluminum foil with a strong, medical-grade acrylic adhesive. Compared to other aluminum foils, AlumaSeal II has less tendency to roll back on itself when removing the backing paper.

- Certified DNase- & RNase- free
- Excellent chemical resistance (e.g. DMSO)
- Easily pierceable with single or multichannel pipettors and robotic probes
- 1.4mil conformable aluminum foil allows visual verification of the seal's security around each well



"For a thicker 2.0mil aluminum foil to prevent tearing upon removal, check out FoilSeal"

