

FoilSeal™ RollMate™

Pre-scored foil in an easy to use roll format

Application:

PCR, light sensitive assays, short-term storage and robotics


Temperature Range:

-20°C to +120°C

CATALOG #	DESCRIPTION	QUANTITY
F2-RL-500	FoilSeal RollMate	500

Pre-scored film on a roll is a more economical and efficient alternative to rapidly seal microplates; faster than applying single sheeted film. Self-standing, non-slip supports hold the pre-scored film roll securely in place for quick and easy bench-top dispensement. The roll unwind tension can be easily adjusted to user preference. Continuous backing liner can be torn off or fed directly into a wastebasket for disposal.

- Certified DNase- & RNase- free
- Pierceable with a pipet tip or robotic probe for direct sample recovery
- 2.0mil aluminum foil removes easily without tearing
- The conformable nature of foil allows for the visual verification of the seal's security around each well
- Not recommended for long-term cold storage



"For an equivalent product in a sheeted format, check out FoilSeal"

